

# Applications

Note 175 | March 2010

## Automated PCR clean-up in 96-well plate and 8-well strip format using the MACHEREY-NAGEL NucleoSpin® 8/96 Extract II kits on the epMotion® 5075 from Eppendorf

Henning Risch, Thomas Zinn, MACHEREY-NAGEL GmbH & Co. KG, Düren, Germany  
Daniel Wehrhahn, Eppendorf AG, Hamburg, Germany

### Abstract

In this application note we demonstrate the integration of the MACHEREY-NAGEL NucleoSpin 8/96 Extract II PCR clean-up kits into the epMotion 5075 VAC automated pipetting system. Depending on the requirements in sample throughput DNA can be purified in 8-well strip format or in a 96-well plate format. The procedure for clean-up of PCR fragments requires minimal set-up time. Fragments from 65 bp up to 10 kbp can be purified. The purified PCR fragments were obtained with high recovery rates and were subjected to downstream applications like sequencing. For 96 samples a processing time of approximately 60 min was achieved.

### Introduction

Purification of amplified PCR fragments is often necessary for several downstream applications like sequencing or microarray spotting. The NucleoSpin Extract II kits allow clean-up of PCR fragments using a silica membrane-based bind-wash-elute method. With the purification procedure undesired components (such as detergents, primers, or PCR additives) are removed efficiently [1].

The MACHEREY-NAGEL NucleoSpin Extract II kits allow for PCR clean-up in either medium throughput 8-well strip format or in the 96-well plate format for higher throughput. Both kits, NucleoSpin 8 Extract II or NucleoSpin 96 Extract II can be used fully automated on the epMotion 5075 automated pipetting system. Amplified PCR products are mixed with a binding buffer and loaded onto the NucleoSpin Extract Binding strips or plate. The PCR fragments are bound reversibly to the silica membrane in a subsequent vacuum binding step. Following washing steps and an ethanol evaporation step the purified PCR fragments can be eluted in water or low salt elution buffer. The purified DNA is suitable for subsequent downstream applications, e.g., sequencing or microarray spotting.

### Materials and Methods

Eppendorf epMotion 5075 VAC  
Vac Frame 2  
Vac Holder  
Reservoir 400 mL  
Collection Plate Adapter for MN Tube Strips (for NucleoSpin 8 Extract II kit only)  
Channeling Plate  
Reservoir Rack with Reagent Reservoirs  
MACHEREY-NAGEL NucleoSpin 8 Extract II kit  
MACHEREY-NAGEL NucleoSpin 96 Extract II kit  
PCR fragments

### Product use limitations and safety information

Please read the MACHEREY-NAGEL NucleoSpin 8 Extract II or NucleoSpin 96 Extract II manuals before performing the method for the first time.

### PCR fragments

PCR products of different sizes (65 bp - 1484 bp) were amplified from a pGEM plasmid in a standard PCR reaction (2 ng template DNA, 0.4 µM forward primer, 0.4 µM reverse primer, 1.5 mM MgCl<sub>2</sub>, 0.2 mM dNTP mix, 2.0 U *Taq* polymerase).

**Determination of yield and purity:**

Purified PCR fragments were dried and analyzed by agarose gel electrophoresis in comparison to not purified PCR reactions. DNA was quantified densitometrically using the Scanalytics ONE D Scan software.

**Agarose gel electrophoresis:**

PCR fragments were analyzed by TAE agarose gel electrophoresis (1 % (w/v) agarose, ethidium bromide stain).

**DNA sequencing**

Purified DNA samples were sequenced using ABI BigDye 3.1 chemistry. Sequence was determined using the ABI 3730 XL DNA Sequencer (MWG sequencing service).

**Table 1:** epMotion 5075 VAC worktable details for NucleoSpin 8 Extract II protocol

| Position | Labware   | Comment  |
|----------|---|--|
| A2       | MN Tube Strips (MN_TP_1200_48)  | elution tubes* (***)   |
| A3       | epT.I.P.S Motion 1000 µL  |  |
| A4       | PCR sample plate (EP_TT_PCR_150)  | 150 µL twin.tec PCR plate  |
| B2       | epT.I.P.S Motion 300 µL   |  |
| B3       | Reagent Reservoirs<br>Position 1: Buffer NT<br>Position 2: Buffer NT3<br>Position 3: empty<br>Position 4: Buffer NE                 | 100 mL reservoir<br>100 mL reservoir<br>empty<br>30 mL reservoir   |
| Vacuum   | NucleoSpin Extract II Binding Strips (MN_FP_8_1400) n Column Holder A**<br>Vacuum Frame 2<br>Reservoir 400 mL with channeling plate | DNA binding strips<br>collar for vacuum manifold<br>collects waste |
| C4       | Vacuum Frame Holder   | Height adapter for vacuum Frame 2                                  |
| T0       | Gripper   |  |
| T1       | TM 1000-8   | 1000 µL 8-channel pipetting tool                                   |
| T2       | TM 300-8  | 300 µL 8-channel pipetting tool                                    |


\*) require Collection Plate Adapter for MN tube strips, see ordering information

\*\*) 8-well strips are inserted into MACHEREY-NAGEL Column Holder A which is a part of the Starter Set A, see ordering information


\*\*\*) 96 well MTP can be used optionally

**Table 2:** epMotion 5075 VAC worktable details for NucleoSpin 96 Extract II protocol

| Position | Labware  | Comment   |
|----------|--|---|
| A2       | Elution Plate (MN_MTP_320)   | MTP for elution   |
| A3       | epT.I.P.S Motion 1000 µL   |   |
| A4       | PCR sample plate (EP_TT_PCR_150)   | 150 µL twin.tec PCR plate   |
| B2       | epT.I.P.S Motion 300 µL  |   |
| B3       | Reagent Reservoirs<br>Position 1: Buffer NT<br>Position 2: Buffer NT3<br>Position 3: Buffer NT3<br>Position 4: Buffer NE | 100 mL reservoir<br>100 mL reservoir<br>100 mL reservoir<br>30 mL reservoir |
| Vacuum   | NucleoSpin Extract II Binding Plate (MN_FP_96_1500)<br>Vacuum Frame 2<br>Reservoir 400 mL with channeling plate          | DNA binding plate<br>collar for vacuum manifold<br>collects waste           |
| C4       | Vacuum Frame Holder  | Height adapter for vacuum Frame 2   |
| T0       | Gripper  |   |
| T1       | TM 1000-8  | 1000 µL 8-channel pipetting tool  |
| T2       | TM 300-8   | 300 µL 8-channel pipetting tool   |



**Fig. 2:** Screenshot from the epMotion Editor showing the setup of the epMotion 5075 VAC worktable for use with the NucleoSpin 96 Extract II




**Fig. 1:** Screenshot from the epMotion Editor showing the setup of the epMotion 5075 VAC worktable for use with the NucleoSpin 8 Extract II kit.


**Result**

**Recovery of PCR fragments**

Total yield of PCR fragments after purification usually depends on the PCR reaction itself, size of the amplified PCR fragment, PCR reaction buffer and PCR reaction additives (e.g. betaine, DMSO, detergents). The NucleoSpin Extract II kits show an excellent performance with a broad range of PCR reaction buffer systems. The kits can be used for purification of PCR fragments down to 65 bp. Figure 3 shows a typical result for the purification of PCR fragments of different sizes. The results for the recovery of the individual PCR fragments are summarized in figure 4.





**Fig. 3 :** Purification of PCR fragments of different sizes. PCR fragments of 65 bp, 164 bp, 252 bp, 359 bp, 459 bp, 645 bp, 782 bp, and 982 bp were purified using the NucleoSpin 96 Extract II kit. Samples were analyzed before purification and following purification on a 1 % TAE agarose gel. The figure shows the purification of samples with the indicated sizes in triplicate. High recovery rates were obtained even with small fragments.



**Fig. 4:** Recovery rates of PCR fragments. PCR fragments with the indicated sizes were purified with the NucleoSpin 96 Extract II kit. Recovery rate was determined densitometrically. Each value represents the average of six individual purifications with the input set to 100 %. High recovery rates of 75 % to 95 % were achieved.

### Sequencing of purified PCR fragments

Primer removal is a critical step for good quality sequencing results of PCR fragments. Incomplete primer removal causes double signals at the beginning of the sequence and may lead to sequencing failure. Using the NucleoSpin 96 Extract II kit allows for purification of small size PCR fragments together with complete removal of PCR primers or primer-dimers. Figure 5 shows the result of a sequencing reaction of a 982 bp PCR fragments. High quality sequencing data together with long reading lengths were obtained.



**Fig. 5:** Automated sequencing of a 982 bp PCR fragment purified with NucleoSpin 96 Extract II kit

### Conclusion

The integration of the MACHEREY-NAGEL NucleoSpin 8 Extract II and 96 Extract II kits into the epMotion 5075 VAC resulted in a flexible system for automated purification of PCR fragments. The system can be used either for low to medium throughput using the 8-well strip based NucleoSpin 8 Extract II kit or for higher throughput using the 96-well based NucleoSpin 96 Extract II kit. Both kits can be used with the same hardware allowing the user to switch between the two methods according to the requirements in sample throughput. PCR fragments of different sizes from 65 bp to 1482 bp were purified with recovery rates of 75 % to 95 %. The purified PCR products are of excellent quality and suitable for downstream applications such as DNA sequencing. Combining the NucleoSpin technology and the epMotion 5075 VAC automated pipetting system forms an attractive and versatile system for the automated clean-up of PCR fragments.

## References

**Eppendorf**

Operating Manual for epMotion 5075

- [1] Birnboim, H.C. & Doly, J. (1979)
- Nucleic Acids Res.*
- 7, 1513-1523

**Macherey-Nagel**

NucleoSpin 8 Extract II kit user manual

NucleoSpin 96 Extract II kit user manual

**Ordering Information Eppendorf**

| Product   | Order no. International | Order no. North America |
|---|-------------------------|-------------------------|
| epMotion® 5075 VAC 230 V (vacuum chamber included)            | 5075 000.164            | n/a                     |
| epMotion® 5075 VAC 120 V (vacuum chamber included)            | n/a                     | 960020014               |
| Collection Plate Adapter MN                                   | 5075 785.064            | 960002571               |
| Channeling Plate  | 5075 794.004            | 960002540               |
| Vac Frame 2   | 5075 785.005            | 960002261               |
| Dispensing tool TM 1000-8                                     | 5280 000.258            | 960001061               |
| Reservoir Rack  | 5075 754.002            | 960002148               |
| Reservoirs 100 mL (10 x 5 reservoirs in bags/case, PCR clean) | 0030 126.513            | 960051017               |
| Reservoirs 30 mL (10 x 5 reservoirs in bags/case, PCR clean)  | 0030 003.993            | 960050100               |

**Ordering Information MACHEREY-NAGEL**

| Product  | Order no.  |
|--|------------|
| NucleoSpin® 8 Extract II kit 12 x 8 preps                  | 740 668    |
| NucleoSpin® 8 Extract II kit 60 x 8 preps                  | 740 668.5  |
| NucleoSpin® 96 Extract II kit 1 x 96 preps                 | 740 658.1  |
| NucleoSpin® 96 Extract II kit 2 x 96 preps                 | 740 658.2  |
| NucleoSpin® 96 Extract II kit 4 x 96 preps                 | 740 658.4  |
| NucleoSpin® 96 Extract II kit 24 x 96 preps                | 740 658.24 |
| Starter Set A (for NucleoSpin 8 Extract II kit only) 1 set | 740 682    |

NucleoSpin® is a trademark of MACHERY-NAGEL GmbH &amp; Co. KG, Dueren, Germany

**MACHEREY-NAGEL**

MACHEREY-NAGEL GmbH & Co. KG  
 Neumann-Neander-Str. 6-8 · D-52355 Dueren, Germany  
 E-mail: tech-bio@mn-net.com · www.mn-net.com

**eppendorf**  
*In touch with life*

Your local distributor: [www.eppendorf.com/worldwide](http://www.eppendorf.com/worldwide)Eppendorf AG • 22331 Hamburg • Germany • Tel. +49 40 53801-0 • Fax +49 40 538 01-556 • E-mail: [eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com)

Eppendorf North America Inc. • 102 Motor Parkway, Suite 410 • Hauppauge, N.Y. 11788-5178 • USA

Tel. +1 516 334 7500 • Toll free phone +1 800 645 3050 • Fax +1 516 334 7506 • E-mail: [info@eppendorf.com](mailto:info@eppendorf.com)**Application Support**Europe, International: Tel: +49 1803 666 789 · E-mail: [support@eppendorf.com](mailto:support@eppendorf.com)North America: Tel: +1 800 645 3050 menu option 2 · E-mail: [techserv@eppendorf.com](mailto:techserv@eppendorf.com)Asia Pacific: Tel: +60 3 8023 6869 · E-mail: [support\\_asiapacific@eppendorf.com](mailto:support_asiapacific@eppendorf.com)