

I. Introduction

This protocol is for use with the PrepX DNA Library Kit, 96 Samples (Cat. No. 640102), which accommodates rapid, walkaway automation of NGS library preparation on the Apollo™ Library Prep system. When run with the PrepX ILMN 32i script, this kit can be used to process **three batches of up to 32 samples per batch** into 220-, 320-, 520-, or 870-base-pair-read DNA libraries suitable for sequencing on Illumina® platforms. **Read this Protocol-At-A-Glance in its entirety before you begin with particular attention paid to the Apollo System Best Practices.**

II. Workflow Overview



Figure 1. Library preparation workflow overview for one batch of up to 32 samples on the Apollo system with the PrepX ILMN DNA Library Kit - B. Blue and purple boxes indicate steps performed on and off the Apollo system, respectively. The italicized text indicates that the post-PCR cleanup requires a separate protocol. Run = the run time in minutes on the Apollo system, if applicable. Total = total time in minutes spent, including thawing of reagents, reagent and equipment setup, heating and cooling of thermal blocks, incubation of reactions, and automated liquid-handling processes, if applicable.

III. List of Components

PrepX DNA Library Kit, 96 Samples (Cat. No. 640102. Also sold separately.)

Box 1. (Store at -20°C)	
PrepX DNA 32i Ligase Buffer, 96	3 x 450 µl
PrepX PCR Master Mix	4 x 700 µl
PrepX Complete 32i Ligase Enzyme	3 x 45 µl
PrepX Intermediate 32i Enzymes Strips Complete, Orange Strips	24 x 4-tube yellow strips
Box 2. (Store at 4°C)	
PrepX Molecular Grade Water	2 x 100 ml
PrepX Cleanup Beads	4 x 6.8 ml
PrepX 2.5M NaCl	10 ml

Additional Materials Required

Sequencing Indexes

- PrepX Complete ILMN Barcodes 1-96 (Cat. No. 640107)

Apollo Consumables

The following consumables must be purchased separately from Takara Bio USA, Inc. and were used to validate protocols and scripts. **Do not make any substitutions.**

Apollo consumables	Cat. No.	Quantity	Usage/32-rxn run
Apollo Piercing Tips	640085	Box of 1,000 tips	8 tips
Apollo Filter Tips	640084	Box of 960 tips	368 tips
Apollo Reservoirs	640087	Box of 100 reservoirs	2 reservoirs
Apollo Microtiter Plates	640083	Box of 25 plates	1 plate
Apollo 1.1 mL MiniTubes	640088	Box of 960 minitubes	96 minitubes
Apollo 0.2 ml PCR 8-Tube Strips, Clear	640082	Box of 125 strips	27 strips
Apollo Caps for 0.2 ml PCR 8-Tube Strips, Clear	640086	Box of 125 strips	27 strips

General lab equipment, reagents, and consumables


- Single-channel pipettes: 10 µl, 20 µl, 200 µl, and 1,000 µl
- Eight-channel pipettes (recommended): 20 µl and 200 µl
- Filter pipette tips: 2 µl, 20 µl, 200 µl, and 1,000 µl
- PCR thermal cycler
- Covaris instrument and related materials for DNA shearing
- DNA LoBind Tubes (Eppendorf, Cat. No. 0030108051)
- 100% ethanol (EtOH; molecular biology grade)

IV. Apollo System Best Practices

- **Read this Protocol-At-A-Glance in its entirety before you begin.**
- **Before beginning library preparation, follow the Illumina index pooling guidelines to select the right barcode adapters for your set of samples.** (Barcode adapter sequences are in Appendix B.) To avoid index read failure due to registration failure, confirm that the combination of barcode adapter sequences you intend to use will maintain color balance for each base of the index read being sequenced. For more information, visit: <https://www.illumina.com/informatics/sample-experiment-management/sequencing-experiment-setup.html>
- Clean the work surfaces, including the retention plates, with 70% ethanol at least once a week.
- **Restart the instrument before every run. Also, between each subprotocol, perform a power cycle by turning the instrument off, waiting 1 min, and then turning it back on.**
- Discard any deformed plastics.
- Separate partial tube strips with scissors and remove resulting plastic overhangs.
- Spin down reagents before placing them on the deck to avoid air bubbles. **Bubbles at the bottoms of tubes must be removed to ensure accurate volume delivery.**
- Ensure plastics are properly seated on the deck surface **with caps/lids removed.** Be sure to push any tubes down completely and evenly prior to installing the metal retention plates.
- Empty the waste box before every run. **An accumulation of tips in the waste box may cause the run to fail.**

V. Protocols

A. Protocol: Sample and Reagent Prep



For each protocol, the corresponding step in the workflow diagram is indicated in with a green outline.

Materials Required

Reagents	Storage conditions	Source
Sample (DNA)	–20°C	User
Plate containing barcode adapters	–20°C	Takara Bio
PrepX DNA 32i Ligase Buffer, 96	–20°C	Takara Bio
PrepX Complete 32i Ligase Enzyme	–20°C	Takara Bio
PrepX Molecular Grade Water	4°C	Takara Bio
PrepX Cleanup Beads	4°C	Takara Bio
PrepX 2.5M NaCl	4°C	Takara Bio
70% ethanol (prepared fresh)	Room temperature	User

NOTE: PrepX Cleanup Beads need to come to room temperature before the container is opened. Therefore, **we strongly recommend preparing ~1-ml aliquots upon receipt** and then refrigerating the aliquots. Individual tubes can be removed for each experiment, allowing them to come to room temperature more quickly (≥ 30 minutes). This will also decrease the chance of bead contamination. Mix well to disperse the beads before adding them to your reactions. The beads are viscous, so pipette slowly.

Apollo consumables	Cat. No.	Quantity	Usage/32-rxn run
Apollo 0.2 ml PCR 8-Tube Strips, Clear	640082	Box of 125 strips	17 strips
Apollo Caps for 0.2 ml PCR 8-Tube Strips, Clear	640086	Box of 125 strips	8 strips
Apollo Microtiter Plates	640083	Box of 25 plates	1 plate

1. Prepare up to 32 strips of 4-well strip tubes by cutting up to 16 8-tube strips in half and trimming the resulting plastic overhangs. The samples and reagents must be distributed in sets of four to be in the correct orientation to process up to 32 samples. **Do not use 8-tube strips instead of 4-tube strips.**
2. For each sample, prepare 1–100 ng of appropriately-sized DNA as follows:
 - a) Quantify DNA using a Qubit Fluorometer or your preferred method.
 - b) Assess the size using a fragment analyzer.
 - c) (Optional) Further fragment the samples using a Covaris instrument according to the parameters in Appendix A and per the manufacturer’s instructions.
 - d) Quantify the DNA again, then prepare the samples at the desired input amount. We recommend using an input amount of 1–100 ng fragmented DNA in a volume of 15 μ l (0.067–6.67 ng/ μ l).

3. Prepare the barcode adapters as follows:
- a) Place the barcode adapter plate in a thermal cycler set to the following annealing conditions:

	Temperature	Time
Step 1	95°C	5 min
Step 2	70°C	15 min

NOTE: This protocol is designed for the generation of up to 8 libraries at a time, so no more than 8 barcodes are used during a given run. Store the plate with the remaining barcode adapters at –20°C until the next use. Do not repeat the annealing process. Move the plate from the thermal cycler to a benchtop and allow it to cool to room temperature.

- b) Add 6 µl of PrepX Molecular Grade Water to each well containing a barcode adapter (96 wells total; use the layout shown in Appendix B, and **always change tips to avoid cross contamination**).
- c) Transfer one batch of barcode adapters into fresh 8-tube strips. Cap and place on ice until used.
4. Prepare the barcode adapter master mix (“Adapter MM”) for **each** of the desired barcode adapters (up to 32) as described in the table below. Add 15 µl of Adapter MM to each tube of a new Apollo 0.2 ml PCR 8-tube strip. Mix thoroughly with a pipette. Cap and keep on ice.

Component	Volume/rxn
Barcode adapter (e.g., well A1)	2 µl
PrepX Molecular Grade Water	13 µl
Total volume	15 µl

5. Prepare the ligation master mix (“Ligation MM”) in one PCR tube as described in the table below. Mix thoroughly with a pipette. Keep on ice until aliquoted.

Component	Volume/rxn	Volume/32 rxns + 1 rxn excess
PrepX DNA 32i Ligase Buffer, 96	12 µl	396 µl
PrepX Complete 32i Ligase Enzyme	1 µl	33 µl
PrepX Molecular Grade Water	2 µl	66 µl
Total volume	15 µl	495 µl

6. On the benchtop, aliquot the reagents into the consumables as described in the table below:

Component	Consumable	Volume/rxn
Sample	New Apollo 8-tube strips	15 µl
Adapter MM	Apollo 8-tube strips (from Step 4)	15 µl
Ligation MM	1 new Apollo 8-tube strip	15 µl
PrepX Molecular Grade Water	1 Apollo reservoir	15 ml
PrepX Cleanup Beads*	New Apollo 8-tube strips	200 µl
PrepX 2.5M NaCl	Apollo microtiter plate, Row 1 (Figure 2)	100 µl
70% ethanol	1 Apollo reservoir	15 ml

*The PrepX Cleanup Beads should be warmed to room temperature and mixed well prior to use.

NOTES:

- The Apollo system is specifically calibrated for the consumables indicated in the table above. Using alternative consumables may cause the run to fail.
- We recommend moving to the Library Synthesis protocol (Section V.B) immediately.

B. Protocol: Library Synthesis



Materials Required

Reagents	Current temperature	Source
Sample	On ice	User
Adapter MM	On ice	Section V.A.
Ligation MM	On ice	Section V.A.
PrepX Molecular Grade Water	Room temperature	Takara Bio
PrepX Cleanup Beads	Room temperature	Takara Bio
PrepX 2.5M NaCl	Room temperature	Takara Bio
70% ethanol	Room temperature	User
PrepX Intermediate 32i Enzymes Strips Complete, Orange Strips	On ice	Takara Bio

Apollo consumables	Cat. No.	Quantity	Usage/32-rxn run
Apollo Piercing Tips	640085	Box of 1,000 tips	8 tips
Apollo Filter Tips	640084	Box of 960 tips	240 tips
Apollo Reservoirs	640087	Box of 100 reservoirs	2 reservoirs
Apollo Microtiter Plates	640083	Box of 25 plates	1 plate
Apollo 1.1 mL MiniTubes	640088	Box of 960 minitubes	96 minitubes
Apollo 0.2 ml PCR 8-Tube Strips, Clear	640082	Box of 125 strips	4 strips
Apollo Caps for 0.2 ml PCR 8-Tube Strips, Clear	640086	Box of 125 strips	4 strips

1. Turn on the instrument or, if the instrument is already on, perform a power cycle by turning the instrument off, waiting 1 min, and then turning it back on.
2. Load Apollo consumables onto the Apollo system work surface according to the layout in Figure 2. First, load the Apollo consumables that do not initially hold reagents (table above). Just before the run, load the Apollo consumables containing reagents, but not samples, onto the system. For information on the deck layout when processing fewer than 32 samples, please see Appendix C.
3. To access the PrepX ILMN 32i script, press **Library Prep > DNA > ILM**. Then, under the PrepX ILM 32i heading, select the desired library insert size (**220 bp**, **320 bp**, **520 bp**, or **870 bp**). The **Cooling** indicator will appear.
4. When the **Cooling** indicator has disappeared, and the **Run** button has appeared, load the samples, reagents, and remaining Apollo consumables onto the Apollo deck following the layout shown in Figure 2 and the instructions shown on the touch screen.

NOTE: Ensure plastics are properly seated on the deck surface with caps carefully removed. Be sure to push any tubes down completely and evenly prior to installing the metal retention plates.

PrepX™ DNA Library Kit, 96 Samples Protocol-At-A-Glance

5. Install the metal retention plates on Blocks 3 and 4.
6. Empty the waste box and remove any used consumables from the system.

NOTE: An accumulation of tips in the waste box may cause the run to fail.

7. Close the instrument door and press **Run**.

NOTE: The run time is 3 hours and 20 minutes.

8. When the run is complete, remove the libraries from Block 3, Row 5 (Products), cap the tubes, and place them on ice. The final DNA library product volume should be ~15 µl per tube.
9. Turn off the instrument.

SAFE STOPPING POINT: If you do not plan to proceed immediately to the Library Amplification protocol, the products in Block 3 can be capped and stored at –20°C for up to one week.

Library Synthesis (32 samples)





Figure 2. Deck layout for library synthesis from 32 samples. For this protocol, use script PrepX_ILMN_32.i. Run time is 3 hours and 20 minutes.

C. Protocol: Library Amplification



Materials Required

Reagents	Current temperature	Source
Products (from Block 3, Row 5)	On ice	Section V.B.
PrepX PCR Master Mix	On ice	Takara Bio
PrepX PCR Primers	On ice	Takara Bio
PrepX Molecular Grade Water	On ice	Takara Bio

Apollo consumables	Cat. No.	Quantity	Usage/32-rxn run
Apollo 0.2 ml PCR 8-Tube Strips, Clear	640082	Box of 125 strips	4 strips
Apollo Caps for 0.2 ml PCR 8-Tube Strips, Clear	640086	Box of 125 strips	4 strips

1. On ice, prepare the Amplification MM in a 1.5-ml tube, per the table below. Mix the Amplification MM by vortexing gently, and then spin down. Aliquot 27 μ l of this master mix into the wells indicated in Figure 2.

Component	Volume/rxn	Volume/32 rxns + 1 rxn excess
PrepX PCR Master Mix	25 μ l	825 μ l
PrepX PCR Primers	2 μ l	66 μ l
Total volume	27 μ l	891 μ l

NOTE: Do not change the volumes of the PrepX PCR Master Mix or the PrepX PCR primers.

2. Use the Amplification MM from Step 1 to prepare the PCR mixture as follows:

Component	Volume/rxn
Amplification Master Mix (from Step 1)	27 μ l
Products (from Block 3, Row 5)	5 ng
PrepX Molecular Grade Water	Varies
Total volume	50 μ l

3. Perform PCR using the following program:

98°C	30 sec
5–10 cycles:	
98°C	10 sec
60°C	30 sec
72°C	30 sec
72°C	300 sec
4°C	forever

NOTE: Depending on the sample type, the PCR conditions may need to be optimized.

4. Spin the tubes to collect the contents, and then bring each reaction volume up to 50 μ l with PrepX Molecular Grade Water. Proceed to the **PrepX PCR Cleanup 32** protocol immediately **OR** PCR products may be stored at 4°C overnight.

Appendix A. Fragmentation Guidelines

Fragment the intact DNA according to the recommended settings below (for a Covaris S220 instrument). Depending on the instrument and sample type used, these settings may need to be optimized. Follow the manufacturer’s instructions if using a different instrument.

Table I. Recommended Covaris S220 Instrument Settings for Desired Library Insert Sizes

Parameter	Recommended settings for library insert size		
	220 bp & 320 bp	520 bp	870 bp
Duty cycle	10%	5%	5%
Intensity	4	3	3
Cycles/burst	200	200	200
Time/cycle	9 sec	10 sec	10 sec
Number of cycles	8	5	4
Total process time	72 sec	50 sec	40 sec

Appendix B. PrepX Complete ILMN Barcodes 1-96

PrepX Complete ILMN Barcodes 1-96 (Cat. No. 640107) allows the preparation of up to 96 Illumina-ready DNA libraries using the PrepX ILMN DNA Library Kit - B. Before beginning library preparation, follow the Illumina index pooling guidelines to select the right barcode adapters for your set of samples. To avoid index read failure due to registration failure, confirm that the combination of barcode adapter sequences you intend to use will maintain color balance for each base of the index read being sequenced. Barcode adapter sequences are listed in Table II. For more information, visit: <https://www.illumina.com/informatics/sample-experiment-management/sequencing-experiment-setup.html>




Figure 3. Contents of PrepX Complete ILMN Barcodes 1-96. The green circle in the diagram represents the PrepX PCR Primers supplied in a tube with a green cap.

PrepX™ DNA Library Kit, 96 Samples Protocol-At-A-Glance

Table II. Barcode adapter sequences for the PrepX Complete ILMN Barcodes 1-96. Barcode sequences are highlighted in red.

Index	Barcode adapter sequence (5'–3')
1	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CTCACG ATCTCGTATGCCGTCTTCTGCTTG
2	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CAGCTT ATCTCGTATGCCGTCTTCTGCTTG
3	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CGCTAC ATCTCGTATGCCGTCTTCTGCTTG
4	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC ATTGTA ATCTCGTATGCCGTCTTCTGCTTG
5	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TCCTAG ATCTCGTATGCCGTCTTCTGCTTG
6	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC AGATCC ATCTCGTATGCCGTCTTCTGCTTG
7	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TTGTCA ATCTCGTATGCCGTCTTCTGCTTG
8	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC AGCTTT ATCTCGTATGCCGTCTTCTGCTTG
9	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC ATCCGC ATCTCGTATGCCGTCTTCTGCTTG
10	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CTGAAA ATCTCGTATGCCGTCTTCTGCTTG
11	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CTGGCC ATCTCGTATGCCGTCTTCTGCTTG
12	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TAATGT ATCTCGTATGCCGTCTTCTGCTTG
13	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC ATTTCT ATCTCGTATGCCGTCTTCTGCTTG
14	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TGTACG ATCTCGTATGCCGTCTTCTGCTTG
15	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CAGTGT ATCTCGTATGCCGTCTTCTGCTTG
16	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GACTCA ATCTCGTATGCCGTCTTCTGCTTG
17	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC ATGACT ATCTCGTATGCCGTCTTCTGCTTG
18	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GTAGGC ATCTCGTATGCCGTCTTCTGCTTG
19	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC AGACCA ATCTCGTATGCCGTCTTCTGCTTG
20	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TCAGCC ATCTCGTATGCCGTCTTCTGCTTG
21	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC ACGGTC ATCTCGTATGCCGTCTTCTGCTTG
22	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TTCGAA ATCTCGTATGCCGTCTTCTGCTTG
23	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC AGGTAC ATCTCGTATGCCGTCTTCTGCTTG
24	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TATCAG ATCTCGTATGCCGTCTTCTGCTTG
25	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CCATGT ATCTCGTATGCCGTCTTCTGCTTG
26	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC ATCGGC ATCTCGTATGCCGTCTTCTGCTTG
27	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TTAGCT ATCTCGTATGCCGTCTTCTGCTTG
28	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GCCATA ATCTCGTATGCCGTCTTCTGCTTG
29	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC AGTGCC ATCTCGTATGCCGTCTTCTGCTTG
30	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CTTGAC ATCTCGTATGCCGTCTTCTGCTTG
31	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CATTAG ATCTCGTATGCCGTCTTCTGCTTG
32	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TCGGAT ATCTCGTATGCCGTCTTCTGCTTG
33	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CGATGT ATCTCGTATGCCGTCTTCTGCTTG
34	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TGACCA ATCTCGTATGCCGTCTTCTGCTTG
35	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GCCAAT ATCTCGTATGCCGTCTTCTGCTTG
36	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CTTGTA ATCTCGTATGCCGTCTTCTGCTTG

PrepX™ DNA Library Kit, 96 Samples Protocol-At-A-Glance

37 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**ATCACG**ATCTCGTATGCCGTCTTCTGCTTG
38 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**TTAGGC**ATCTCGTATGCCGTCTTCTGCTTG
39 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**ACTTGA**ATCTCGTATGCCGTCTTCTGCTTG
40 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**GATCAG**ATCTCGTATGCCGTCTTCTGCTTG
41 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**TAGCTT**ATCTCGTATGCCGTCTTCTGCTTG
42 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**GGCTAC**ATCTCGTATGCCGTCTTCTGCTTG
43 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**AGTCAA**ATCTCGTATGCCGTCTTCTGCTTG
44 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**AGTTCC**ATCTCGTATGCCGTCTTCTGCTTG
45 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**ATGTCA**ATCTCGTATGCCGTCTTCTGCTTG
46 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**CCGTCC**ATCTCGTATGCCGTCTTCTGCTTG
47 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**GTAGAG**ATCTCGTATGCCGTCTTCTGCTTG
48 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**GTCCGC**ATCTCGTATGCCGTCTTCTGCTTG
49 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**GTGAAA**ATCTCGTATGCCGTCTTCTGCTTG
50 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**GTGGCC**ATCTCGTATGCCGTCTTCTGCTTG
51 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**GT TTCG**ATCTCGTATGCCGTCTTCTGCTTG
52 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**CGTACG**ATCTCGTATGCCGTCTTCTGCTTG
53 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**GAGTGG**ATCTCGTATGCCGTCTTCTGCTTG
54 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**GGTAGC**ATCTCGTATGCCGTCTTCTGCTTG
55 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**ACTGAT**ATCTCGTATGCCGTCTTCTGCTTG
56 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**ATGAGC**ATCTCGTATGCCGTCTTCTGCTTG
57 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**ATTCCT**ATCTCGTATGCCGTCTTCTGCTTG
58 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**CAAAG**ATCTCGTATGCCGTCTTCTGCTTG
59 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**CAACTA**ATCTCGTATGCCGTCTTCTGCTTG
60 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**CACCGG**ATCTCGTATGCCGTCTTCTGCTTG
61 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**CACGAT**ATCTCGTATGCCGTCTTCTGCTTG
62 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**CACTCA**ATCTCGTATGCCGTCTTCTGCTTG
63 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**CAGGCG**ATCTCGTATGCCGTCTTCTGCTTG
64 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**CATGGC**ATCTCGTATGCCGTCTTCTGCTTG
65 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**CATTTT**ATCTCGTATGCCGTCTTCTGCTTG
66 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**CCAACA**ATCTCGTATGCCGTCTTCTGCTTG
67 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**CGGAAT**ATCTCGTATGCCGTCTTCTGCTTG
68 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**CTATAC**ATCTCGTATGCCGTCTTCTGCTTG
69 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**CTCAGA**ATCTCGTATGCCGTCTTCTGCTTG
70 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**GACGAC**ATCTCGTATGCCGTCTTCTGCTTG
71 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**TAATCG**ATCTCGTATGCCGTCTTCTGCTTG
72 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**TACAGC**ATCTCGTATGCCGTCTTCTGCTTG
73 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**TATAAT**ATCTCGTATGCCGTCTTCTGCTTG
74 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**TCATTC**ATCTCGTATGCCGTCTTCTGCTTG
75 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**TCCCGA**ATCTCGTATGCCGTCTTCTGCTTG
76 GATCGGAAGAGCACACGTCTGAACTCCAGTCAC**TCGAAG**ATCTCGTATGCCGTCTTCTGCTTG

PrepX™ DNA Library Kit, 96 Samples Protocol-At-A-Glance

77	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TTGACT ATCTCGTATGCCGTCTTCTGCTTG
78	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GGAACT ATCTCGTATGCCGTCTTCTGCTTG
79	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TGACAT ATCTCGTATGCCGTCTTCTGCTTG
80	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GGACGG ATCTCGTATGCCGTCTTCTGCTTG
81	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CTCTAC ATCTCGTATGCCGTCTTCTGCTTG
82	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GCGGAC ATCTCGTATGCCGTCTTCTGCTTG
83	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TTTCAC ATCTCGTATGCCGTCTTCTGCTTG
84	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GGCCAC ATCTCGTATGCCGTCTTCTGCTTG
85	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CGAAAC ATCTCGTATGCCGTCTTCTGCTTG
86	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC ACACTC ATCTCGTATGCCGTCTTCTGCTTG
87	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GCTACC ATCTCGTATGCCGTCTTCTGCTTG
88	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC ATCAGT ATCTCGTATGCCGTCTTCTGCTTG
89	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC GCTCAT ATCTCGTATGCCGTCTTCTGCTTG
90	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC AGGAAT ATCTCGTATGCCGTCTTCTGCTTG
91	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TTTTTG ATCTCGTATGCCGTCTTCTGCTTG
92	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TAGTTG ATCTCGTATGCCGTCTTCTGCTTG
93	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CCGGTG ATCTCGTATGCCGTCTTCTGCTTG
94	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC ATCGTG ATCTCGTATGCCGTCTTCTGCTTG
95	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC TGAGTG ATCTCGTATGCCGTCTTCTGCTTG
96	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC CGCCTG ATCTCGTATGCCGTCTTCTGCTTG

Appendix C. Deck Setup for Additional Reaction Sizes

Using the PrepX ILMN 8 and PrepX Cleanup 8 scripts, the Apollo system can run up to eight samples per batch. If you are performing fewer than eight reactions, please use a specific deck layout based on the number of samples (Table VI).

Table III. Deck Layout Options for Various Sample Numbers


Setting up less than 8 samples	
# of samples	Columns to load
1	D or E*
2	D, E
3	C, D, E or D, E, F
4	C, D, E, F
5	B, C, D, E, F or C, D, E, F, G
6	B, C, D, E, F, G
7	A, B, C, D, E, F, G or B, C, D, E, F, G, H

*The user can choose between a deck setup centered around column D or E. The scripts are designed to accommodate either layout.


Appendix D. Library Validation

Figures 4 and 5 show example traces from pre-PCR and post-PCR library outputs using the 220-, 320-, 520-, or 870-bp size selection options.


Insert size: 220 bp



Insert size: 320 bp



Insert size: 520 bp



Insert size: 870 bp






Figure 4. Bioanalyzer traces of pre-PCR library outputs. Covaris-fragmented human DNA was fragmented via sonication, processed through the Protocol: Library Amplification (Section V.C) and validated using an Agilent 2100 Bioanalyzer and Agilent’s High Sensitivity DNA Kit. A no-template control is shown in blue.


Insert size: 220 bp



Insert size: 320 bp



Insert size: 520 bp



Insert size: 870 bp




Figure 5. Bioanalyzer traces of post-PCR library outputs. The pre-PCR library outputs shown in Figure 5 were amplified in a commercial thermal cycler for a total of five cycles, then cleaned up using the PrepX PCR Cleanup 32 Protocol-At-A-Glance. No-template controls are shown in blue.

Contact Us	
Customer Service/Ordering	Technical Support
tel: 800.662.2566 (toll-free)	tel: 800.662.2566 (toll-free)
fax: 800.424.1350 (toll-free)	fax: 800.424.1350 (toll-free)
web: takarabio.com	web: takarabio.com
e-mail: ordersUS@takarabio.com	e-mail: technical_support@takarabio.com

Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

©2018 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

This document has been reviewed and approved by the Quality Department.