

ApoAlert® Annexin V User Manual



Clontech

United States/Canada
800.662.2566

Asia Pacific
+1.650.919.7300

Europe
+33.(0)1.3904.6880

Japan
+81.(0)77543.6116

Clontech Laboratories, Inc.
A Takara Bio Company
1290 Terra Bella Ave.
Mountain View, CA 94043
Technical Support (US)
E-mail: tech@clontech.com
www.clontech.com

Cat. Nos. 630109, 630110
630201 & 630202
PT3050-1 (PR983322)
Published August 2009

Table of Contents

I. Introduction	3
II. List of Components	4
III. Additional Materials Required	5
IV. Annexin V Protocol	6
A. Incubation of Cells with Annexin V	6
B. Quantification by Flow Cytometry	6
C. Detection by Fluorescence Microscopy	6
V. References	7

Contact Us For Assistance

Customer Service/Ordering:	Technical Support:
Telephone: 800.662.2566 (toll-free)	Telephone: 800.662.2566 (toll-free)
Fax: 800.424.1350 (toll-free)	Fax: 650.424.1064
Web: www.clontech.com	Web: www.clontech.com
E-mail: orders@clontech.com	E-mail: tech@clontech.com

I. Introduction

Apoptosis, or programmed cell death, plays a fundamental role in many normal biological processes as well as several disease states (Wyllie *et al.*, 1980; Ellis *et al.*, 1991; Cohen *et al.*, 1992; Thompson, 1995; Nicholson, 1996). Apoptosis can be induced by various stimuli that all produce the same end result: systematic and deliberate cell death.

One method for studying apoptosis detects changes in the position of phosphatidylserine (PS) in the cell membrane. In non-apoptotic cells, most PS molecules are localized at the inner layer of the plasma membrane, but soon after inducing apoptosis, PS redistributes to the outer layer of the membrane, and becomes exposed to the extracellular environment (Martin *et al.*, 1995; Fadok *et al.*, 1992a; 1992b; 1993). PS translocation precedes other apoptotic events, thus allowing early detection of apoptosis.

Exposed PS can be easily detected with annexin V, a 35.8-kDa protein that has a strong affinity for PS. The specificity of the annexin V-PS interaction is reflected by competitive binding studies in which annexin V binding to apoptotic cells was specifically inhibited by PS liposomes, but was unaffected by liposomes containing other phospholipids (Martin *et al.*, 1995). The assay is nonenzymatic, and can be used with flow cytometry or fluorescence microscopy.

PS externalization occurs independently of the stimulus used. For example, treatment with antibodies to the Fas cell surface receptor induces apoptosis in certain cell types like cells of the immune system (Yonehara *et al.*, 1989; Rouvier *et al.*, 1993; Griffith, 1995). Fas-induced apoptosis in Jurkat T cells is profiled with the Annexin V Assay in Figure 1.




Figure 1. Annexin V-FITC staining of apoptosis in Jurkat T cells. Apoptosis was induced with 200 ng/ml anti-Fas (IgM, CH11). At the indicated times, cells were incubated for 5 min with 1.0 µg/ml annexin V-FITC in binding buffer, and 5,000 cells were examined by flow cytometry analysis. For details, see Martin *et al.*, 1995. Data reproduced with author's permission.

II. List of Components

Store all components at 4°C.

Cat. No.	Cat. No.	
630109	630110	
50 assays	200 assays	
• 250 µl	1 ml	Annexin V-FITC (20 µg/ml in Tris-NaCl)
• 2 x 25 ml	2 x 100 ml	1X Binding Buffer
• 500 µl	2 ml	Propidium iodide (50 µg/ml in 1X binding buffer)

Cat. No. 630201

- 500 µl **Annexin V-FITC** (100 µg/ml)

Cat. No. 630202

- 100 ml **10X Annexin V Binding Buffer**

Notes:

- Under acid conditions, sodium azide yields hydrazoic acid, which is extremely toxic. Azide compounds should be discarded by flushing with running water. We recommend these precautions to avoid deposits in metal piping, in which explosive condition can develop. If the event of skin or eye contact, wash excessively with water. Discard all samples and material coming in contact with sodium azide with due precaution.
- Do not expose reagent to strong light during storage and incubation.
- Microbial contamination of reagents may cause incorrect results.
- Propidium iodide is toxic. Handle with extreme caution.

III. Additional Materials Required

Annexin V Protocol:

- **PBS**

	<u>Final conc.</u>	<u>To prepare 2 L of solution:</u>
Na_2HPO_4	58 mM	16.5 g
NaH_2PO_4	17 mM	4.1 g
NaCl	68 mM	8.0 g

Dissolve components in 1.8 L of H_2O . Adjust to pH 7.4 with 0.1 N NaOH. Add H_2O to a final volume of 2 L. Store at room temperature.

- **Centrifuge** for collecting cells

IV. Annexin V Protocol

PLEASE READ THE ENTIRE PROTOCOL BEFORE BEGINNING

Because apoptosis is a rapid and dynamic process, we recommend performing analysis immediately after staining.

If you plan to fix your cells, incubate them with annexin V **before** fixation, because cell membrane disruption can allow annexin V to bind to PS on the inner surface of the cell membrane. Rinse unbound annexin V with binding buffer before fixation.

A. Incubation of cells with Annexin V

1. Induce apoptosis by a desired method.
 - a. For nonadherent cells, proceed to Step 2.
 - b. For **flow cytometry** with adherent cells, trypsinize cells. Gently wash the cells once with serum-containing media before incubation with annexin V. Proceed to Step 2.
 - c. For **fluorescence microscopy** with adherent cells, grow cells directly on a glass slide and proceed to Step 2.
2. Rinse 1×10^5 – 1×10^6 cells with 1X Binding Buffer.
3. Resuspend the cells in 200 μ l of 1X Binding Buffer.
4. Add 5 μ l of Annexin V and [optional] 10 μ l of Propidium Iodide.
5. Incubate at room temperature for 5–15 min in the dark.

Proceed to B or C below, depending on analysis method.

B. Quantification by Flow Cytometry

Analyze cells by flow cytometry using a single laser emitting excitation light at 488 nm. Use Binding Buffer to bring the reaction volume to at least 500 μ l for flow cytometry analysis.

Note: The signal generated by Annexin V can be detected in the FITC signal detector, and the signal generated by PI can be monitored by the detector reserved for phycoerythrin emission.

C. Detection by Fluorescence Microscopy

1. For analyzing nonadherent cells, place the cell suspension from Step A.5 on a glass slide. Cover with a glass coverslip.

For analyzing adherent cells, following incubation (A.2–5), invert coverslip on glass slide and visualize cells. After incubating with annexin V, the cells can be washed and fixed in 2% formaldehyde.
2. Observe the cells under a fluorescence microscope using a dual filter set for FITC & rhodamine.

Note: Cells with bound annexin V will show green staining in the plasma membrane. Cells that have lost membrane integrity will show red staining (PI) throughout the cytoplasm and a halo of green staining on the cell surface (plasma membrane).

V. References

- Cohen, J. J., Duke, R. C., Fadok, V. A. & Sellins, K. S. (1992) Apoptosis and cell death in immunity. *Ann. Rev. Immunol.* **10**:267–293.
- Ellis, R. E., Yuan, J.Y. & Horvitz, H. R. (1991) Mechanisms and function of cell death. *Ann. Rev. Cell Biol.* **7**:663–698.
- Fadok, V. A., Voelker, D. R., Campbell, P. A., Cohen, J. J., Bratton, G. L. & Henson P. M. (1992a) Exposure of phosphatidylserine on the surface of apoptotic lymphocytes triggers specific recognition and removal by macrophages. *J. Immunol.* **148**:2207–2216.
- Fadok, V. A., Savill, J. S., Haslett, D. L., Bratton, G. L., Doherty, P. A., Campbell, P. A. & Henson P. M. (1992b) Different populations of macrophages use either the vitronectin receptor or the phosphatidylserine receptor to recognize and remove apoptotic cells. *J. Immunol.* **149**:4029–4035.
- Fadok, V. A., Lazlo, D. J., Noble, P.W., Weinstein, L., Riches, D.W. H. & Henson P. M. (1993) Particle digestibility is required for induction of the phosphatidylserine recognition mechanism used by murine macrophages to phagocytose apoptotic cells. *J. Immunol.* **151**:4274–4285.
- Griffith, T. S. (1995) Fas ligand-induced apoptosis as a mechanism of immune privilege. *Science* **270**:1189–1192.
- Koester, S. K., Roth, P., Mikulka, W. R., Schlossman, S. F., Zhang, C. & Bolton, W. E. (1997) Monitoring early cellular responses in apoptosis is aided by the mitochondrial membrane protein-specific monoclonal antibody APO2.7. *Cytometry* **29**:306–312.
- Koester, S. K., Schlossman, S. F., Zhang, C., Decker S. J. & Bolton, W. E. (1998) APO2.7 defines a shared apoptic-necrotic pathway in a breast tumor hypoxia model. *Cytometry* **33**:324–332.
- Martin, S. J., Reutelingsperger, C. P., McGahon, A. J., Rader, J. A., van Schie, R. C., LaFace, D. M. & Green, D. R. (1995) Early redistribution of plasma membrane phosphatidylserine is a general feature of apoptosis regardless of the initiating stimulus: inhibition by overexpression of Bcl-2 and Abl. *J. Exp. Med.* **182**:1545–1556.
- Nicholson, D. W. (1996) ICE/CED3-like proteases as therapeutic targets for the control of inappropriate apoptosis. *Nature Biotechnol.* **14**:297–301.
- Rouvier, E., Luciani, M.F. & Golstein, P. (1993) Fas involvement in Ca²⁺-independent T cell-mediated cytotoxicity. *J. Exp. Med.* **177**:195–200.
- Seth, A., Zhang, C., Letvin, N. L. & Schlossman, S. F. (1997) Detection of apoptotic cells from peripheral blood of HIV-infected individuals using a novel monoclonal antibody. *Cytometry* **33**:324–332.
- Thompson, C. B. (1995) Apoptosis in the pathogenesis and treatment of disease. *Science* **267**:1456–1462.
- Wyllie, A. H., Kerr, J. F. & Currie, A. R. (1980) Cell death: The significance of apoptosis. *Intl. Rev. Cytol.* **68**:251–306.
- Yonehara, S., Ishii, A. & Yonehara, M. (1989) A cell-killing monoclonal antibody (anti-Fas) to a cell-surface antigen co-downregulated with the receptor of tumor necrosis factor. *J. Exp. Med.* **169**:1747–1756.
- Zhang, C., Ao, Z., Seth A., & Schlossman, S. F. (1996) A mitochondrial membrane protein defined by a novel monoclonal antibody is preferentially detected in apoptotic cells. *J. Immunol.* **157**:3980–3987.

Notes

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

These products are the subjects of pending U.S. and foreign patents.

Not-For-Profit-Entities: Orders may be placed in the normal manner by contacting your local representative or Clontech Customer Service at 800.662.2566. At its discretion, Clontech grants not-for-profit research entities a worldwide, non-exclusive, royalty-free, limited license to use this product for non-commercial life science research use only. Such license specifically excludes the right to sell or otherwise transfer this product or its components to third parties. Any other use of this product will require a license from Clontech. For license information please contact the licensing hot line at 800.662.2566, extension 7816 or by e-mail at licensing@clontech.com.

For-Profit entities that wish to use this product in non-commercial or commercial applications are required to obtain a license from Clontech. For license information please contact the licensing hot line at 650.919.7320, or by e-mail at licensing@clontech.com.

These products are covered by U.S. Patent Nos. 5,491,084 and 6,146,826 issued to Columbia University.

These products are sold under license from Columbia University. Rights to use this product are limited to internal research use only; **NOT FOR DIAGNOSTIC OR THERAPEUTIC USE IN HUMANS OR ANIMALS**. No other rights are conveyed. Inquiry into the availability of a license to a broader rights or the use of this product for commercial purposes (e.g., use in manufacturing processes or in the development, screening or discovery of products for therapeutic, diagnostic or prophylactic purposes) should be directed to Columbia University, Columbia Innovation Enterprises, Science and Technology Ventures, 363 Engineering Terrace, Mail Code 2206, 500 West 120th Street, New York, New York 10027, USA.

Clontech, Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc., unless noted otherwise. Clontech is a Takara Bio Company. ©2009 Clontech Laboratories, Inc.