

LD-Insert Screening Amplimer Sets Protocol-at-a-Glance (PT1579-2)

This Protocol-at-a-Glance is provided for your convenience, but is not intended for first-time users. Please read the User Manual before using this abbreviated protocol.



Clontech

United States/Canada
800.662.2566

Asia Pacific
+1.650.919.7300

Europe
+33.(0)1.3904.6880

Japan
+81.(0)77.543.6116

Clontech Laboratories, Inc.
A Takara Bio Company
1290 Terra Bella Ave.
Mountain View, CA 94043
Technical Support (US)
E-mail: tech@clontech.com
www.clontech.com

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc.
Clontech is a Takara Bio Company. ©2006
(PR9Y350; published 08 December 1999)

cDNA Insert Screening

Uses:

 λ Libraries

(λgt10, λgt11, λACT, & λTriplEx™)

Plasmid Libraries

(pACT, pACT2, pGAD 10, pGAD GH, pGAD GL, pB42 AD, & pTriplEx)

Template Prep:

E. coli: Pick plaque and transfer to 25 μl deionized H₂O. Vortex. Use 5 μl per reaction.*E. coli*: Pick colony and transfer to 25 μl deionized H₂O. Boil samples 5 min. Use 5 μl per reaction.*Yeast*: See User Manual, Section IV.D. Use 1–5 μl per reaction.

PCR Set-up:

1. Thaw all components on ice. Mix each component thoroughly before use.
2. Prepare PCR Master Mix (prepare enough mix for all PCR reactions + 1):

	1 rxn	(n+1)rxn	Expt.
PCR-grade deionized H ₂ O	36 μl	36(n+1) μl	___ μl
10X Advantage® 2 PCR Buffer	5 μl	5(n+1) μl	___ μl
5' LD Amplimer	1 μl	1(n+1) μl	___ μl
3' LD Amplimer	1 μl	1(n+1) μl	___ μl
50X dNTP Mix (10 mM each)	1 μl	1(n+1) μl	___ μl
Advantage 2 Polymerase Mix (50X)*	1 μl	1(n+1) μl	___ μl

* If you are not using Advantage 2 Polymerase Mix, use the recipe above if Mg⁺⁺ is included in the 10X reaction buffer, or use the recipe for the genomic PCR Master Mix if Mg⁺⁺ is not included in the buffer.

3. Mix components by vortexing, and spin briefly to collect contents at bottom of tube.
4. For each PCR reaction, combine 45 μl of PCR Master Mix and 5 μl of template DNA (use 1–5 μl of yeast-isolated plasmid). Use 5 μl of deionized H₂O as a negative control.
5. Spin tubes briefly to collect contents.
6. Add 1–2 drops of mineral oil to each tube to prevent evaporation.

PCR:

7. Commence thermal cycling in a Applied Biosystems DNA Thermal Cycler Model 480. These are general guidelines—optimal parameters may vary with different thermal cyclers, templates, and other experimental variables.

Target Size	Cycle Parameters	Target Size	Cycle Parameters
< 5 kb:	<ul style="list-style-type: none"> • 94°C for 1 min • 25–35 cycles 94°C 30 sec^a 68°C 3 min • 68°C for 3 min^b • Soak at 15°C 	5–9 kb:	<ul style="list-style-type: none"> • 94°C for 1 min • 25–35 cycles 94°C 30 sec^a 68°C 6 min • 68°C for 6 min^b • Soak at 15°C

^a Use shortest possible denaturation time.

^b Optional: may reduce background.

Analysis:

8. Transfer 5 μl of the PCR product to a fresh tube, and add 1 μl of stop/loading buffer.
9. Electrophorese samples on a suitable agarose/EtBr gel. (We recommend 1.5% agarose for inserts of 0.3–1.5 kb, 1.2% agarose for inserts of 0.5–10 kb, and 0.8% agarose for inserts >5 kb.)

Genomic & Cosmid Insert Screening

Uses: λ Libraries(EMBL3 & EMBL3 SP6/T7 [also known as λ GEM[®] II])Cosmid Libraries

(pWE15)

Template Prep:*E. coli*: Pick plaque and transfer to 25 μ l deionized H₂O. Vortex. Use 5 μ l/PCR rxn.*E. coli*: Standard alkaline lysis method. Resuspend DNA to ~8 ng/ μ l (5-fold more dilute than standard miniprep DNA). Use 5 μ l/PCR rxn.**Set-up PCR:**

1. Thaw all components on ice. Mix each component thoroughly before use.
2. Prepare PCR Master Mix (prepare enough mix for all PCR reactions +1):

	1 rxn	(n+1)rxn	Expt.
PCR-grade deionized H ₂ O	33.8 μ l	33.8(n+1) μ l	___ μ l
10X <i>Tth</i> PCR reaction buffer	5 μ l	5(n+1) μ l	___ μ l
25 mM Mg(OAc) ₂	2.2 μ l	2.2(n+1) μ l	___ μ l
5' LD Amplimer	1 μ l	1(n+1) μ l	___ μ l
3' LD Amplimer	1 μ l	1(n+1) μ l	___ μ l
50X dNTP Mix (10 mM each)	1 μ l	1(n+1) μ l	___ μ l
Advantage <i>Tth</i> Polymerase Mix (50X)*	1 μ l	1(n+1) μ l	___ μ l

* If you are not using Advantage *Tth* Polymerase Mix, use the recipe above if Mg⁺⁺ is not included in the 10X reaction buffer, or use the recipe for the cDNA PCR Master Mix if Mg⁺⁺ is included in the buffer.

3. Mix components by vortexing, and spin briefly to collect contents at bottom of tube.
4. For each PCR reaction, combine 45 μ l of PCR Master Mix and 5 μ l of template DNA. Use 5 μ l of deionized H₂O as a negative control.
5. Spin tubes briefly to collect contents.
6. Add 1–2 drops of mineral oil to each tube to prevent evaporation.

PCR:

7. Commence thermal cycling in a Applied Biosystems DNA Thermal Cycler Model 480. These are general guidelines—optimal parameters may vary with different thermal cyclers, templates, and other experimental variables.

<u>Target Size</u>	<u>Cycle Parameters</u>	<u>Target Size</u>	<u>Cycle Parameters</u>
< 5 kb:	<ul style="list-style-type: none"> • 94°C for 1 min • 25–35 cycles <li style="padding-left: 20px;">94°C 30 sec^a <li style="padding-left: 20px;">68°C 3 min • 68°C for 3 min^b • Soak at 15°C 	10–20 kb:	<ul style="list-style-type: none"> • 94°C for 1 min • 25–35 cycles <li style="padding-left: 20px;">94°C 30 sec^a <li style="padding-left: 20px;">68°C 12 min • 68°C for 12 min^b • Soak at 15°C
<u>Target Size</u>	<u>Cycle Parameters</u>	<u>Target Size</u>	<u>Cycle Parameters</u>
5–9 kb:	<ul style="list-style-type: none"> • 94°C for 1 min • 25–35 cycles <li style="padding-left: 20px;">94°C 30 sec^a <li style="padding-left: 20px;">68°C 6 min • 68°C for 6min^b • Soak at 15°C 	20–40 kb:	<ul style="list-style-type: none"> • 94°C for 1 min • 25–35 cycles <li style="padding-left: 20px;">94°C 15 sec^a <li style="padding-left: 20px;">68°C 22 min • 68°C for 22 min^b • Soak at 15°C

^a Use shortest possible denaturation time.

^b Optional: may reduce background.

Analysis:

8. Transfer 5 μ l of the PCR product to a fresh tube and add 1 μ l of stop/loading buffer.
9. Electrophorese samples on a suitable agarose/EtBr gel. (We recommend 1.5% agarose for inserts of 0.3–1.5 kb, 1.2% agarose for inserts of 0.5–10 kb, and 0.8% agarose for inserts >5 kb.)

© 1999, Clontech Laboratories, Inc.