

Terra™ Direct PCR Polymerase Mix Protocol-At-A-Glance

Please read the User Manual before using this Protocol-at-a-Glance. This abbreviated protocol is provided for your convenience, but is not intended for first-time users.

1. Terra Direct PCR Polymerase Mix is designed to amplify targets from whole tissues and tissue extracts. Use the sample amount—and/or preparation method—that is appropriate for your sample type.
2. On ice, combine the reagents shown in Table I in a PCR tube. Make sure you use enough of each reagent for all of your reactions plus an additional reaction to accommodate pipetting errors.

Table 1. Recommended Reagent Volumes

| Reagent | Volume | Final concentration |
|---|----------|---------------------|
| 2X Terra Direct PCR Buffer (with Mg ²⁺ , dNTP) | 25 ml | 1X |
| Primer 1 | 15 pmol | 0.3 μM |
| Primer 2 | 15 pmol | 0.3 μM |
| Tissue Sample/Extract | ≤ 5 μl * | |
| Terra Direct PCR Polymerase Mix | 1 μl | 1.25 U |
| Sterile Water | to 50 μl | |
| Total volume per reaction | 50 μl † | |

* When using: 1) crude extracts from various animal tissues, add ≤ 5 μl; 2) blood treated with EDTA or heparin, add ≤ 5 μl; 3) mouse tail, add ≤ 1 mm; 4) mouse ear, add ≤ 1.5 mm²; or 5) plant leaves (e.g., tomato or spinach), add ≤ 1.2 mm diameter disk.

† For 25 μl reactions, be sure to add only half the amount indicated for each reagent.

3. Mix the contents of each tube by tapping the bottom of the tube, then centrifuge briefly.
4. Program your thermal cycler with the following cycling conditions:

- **3-Step PCR** (for amplification of standard targets < 2 kb):

| | | | |
|------|----------|---|--------------|
| 98°C | 2 min* | } | 30–40 cycles |
| 98°C | 10 sec | | |
| 60°C | 15 sec | | |
| 68°C | 1 min/kb | | |

* The initial denaturation step must be performed at 98°C for 2 min in order to denature the hot start antibody.

- **2-Step PCR** (for amplification of targets that are GC-rich or ≥ 2 kb):

| | | | |
|------|----------|---|--------------|
| 98°C | 2 min* | } | 30–40 cycles |
| 98°C | 10 sec | | |
| 68°C | 1 min/kb | | |

* The initial denaturation step must be performed at 98°C for 2 min in order to denature the hot start antibody.

5. Post-PCR Considerations

- Use TAE running buffer when visualizing your PCR products by agarose gel electrophoresis.
- Add proteinase K to the gel loading buffer.
 - a. Add 5 μl of Proteinase K to 50 μl of 5–6X loading buffer.
 - b. Before loading your samples onto a gel, add 15 μl of the loading buffer-proteinase K mixture to the entire 50 μl PCR reaction, or add 1 μl of the loading buffer-proteinase K mixture to 4 μl of the PCR reaction.

| Contact Us | |
|--|--|
| Customer Service/Ordering | Technical Support |
| tel: 800.662.2566 (toll-free) | tel: 800.662.2566 (toll-free) |
| fax: 800.424.1350 (toll-free) | fax: 800.424.1350 (toll-free) |
| web: takarabio.com/service | web: takarabio.com/support |
| e-mail: ordersUS@takarabio.com | e-mail: technical_support@takarabio.com |

Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements

© 2024 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

This document has been reviewed and approved by the Quality Department.