

I. List of Components

- 2X Terra Direct PCR Red Dye Premix
- dH₂O (PCR-grade)
- 2 tubes Proteinase K (100 µl/tube)

II. Additional Materials Required

- Gene-specific PCR primers (with $T_m \geq 60^\circ\text{C}$; avoid using primers that contain inosine)
- PCR reaction tubes or plate
- Aerosol-resistant pipette tips preferably equipped with hydrophobic filters.

III. General Considerations

The red dye in the premix serves two functions as a visual aid. First, the colored dye aids positioning during the loading process. Second, the dye separates into two colored components to aid in monitoring the progress of electrophoresis. Follow the recommendations in Option A (below) if you are performing direct PCR with whole-tissue samples and Option B if you are performing PCR with crude tissue extracts.

NOTE: While there is a tolerance in the amount of material added, the maximum recommended amounts listed here are more than sufficient to obtain good signal. In fact, *less sample may provide a higher signal.*

Option A

Direct PCR with Whole-Tissue Samples

Use the following as a guide to help you determine the appropriate amount of whole-tissue sample to use in your direct PCR reaction.

When using:

- blood treated with EDTA or heparin, add ≤ 5 µl of blood ***directly to*** your PCR reaction.
- mouse tail biopsies, add ≤ 1 mm of tissue ***directly to*** your PCR reaction.
- mouse ear biopsies, add ≤ 1.5 mm² of tissue ***directly to*** your PCR reaction.
- plant leaves (e.g., tomato or spinach), add a ≤ 1.2 mm diameter disc ***directly to*** your PCR reaction.

Option B

PCR with Crude Tissue Extracts

For tissue extracts, you may use any sample preparation method that is appropriate for your sample type. We have found the following method works well for mouse tissue extracts:

- Add 180 µl of 50 mM NaOH to 5–10 mg of mouse tissue (e.g., tail, liver, spleen, thymus, or brain) and incubate for ten minutes at 95°C.
- Neutralize the extract by adding 20 µl of 1M Tris-HCl (pH 8.0).
- Add ≤ 5 µl of the crude extract to the PCR reaction (see Table 1, below).

IV. Protocol

A. PCR Reaction Set-Up

1. Prepare the PCR reaction on ice by adding each component indicated in Table 1.
2. Briefly spin the tube in a microcentrifuge, and begin thermal cycling using the guidelines provided below.

Table 1. Recommended Reagent Volumes

Reagent	Amount	Final concentration
2X Terra Direct PCR Red Dye Premix	25 µl	1X
Primer 1	15 pmol	0.3 µM
Primer 2	15 pmol	0.3 µM
Tissue Sample/Extract	≤ 5 µl ^a	
dH ₂ O	to 50 µl ^b	
Total volume per reaction	50 µl ^c	

^a See the ‘General Considerations’ section, above, for suggested amounts of different sample types.

^b In the PCR reaction, the **final concentration** of Mg²⁺ is 2 mM and the final concentration of each dNTP is 400 mM.

^c For 25 µl reactions, be sure to add only half the amount indicated for each reagent.

B. Recommended Cycling Conditions

Use the following cycling conditions when setting up your initial experiments. These are general guidelines—the optimal conditions may vary.

3-Step PCR (For amplification of standard targets)				2-Step PCR (For amplification of GC-rich targets)			
98°C	2 min*] 30–40 cycles		98°C	2 min*] 30–40 cycles	
98°C	10 sec			98°C	10 sec		
60°C	15 sec			68°C	1 min/kb		
68°C	1 min/kb						

* The initial denaturation step must be performed at 98°C for 2 min in order to denature the hot start antibody.

C. Post-PCR Considerations

- Use TAE running buffer for agarose gel electrophoresis. TBE is not recommended as it causes spreading of the DNA bands toward the bottom of the gel. Load 3–5 µl of each reaction on the gel after cycling.
- PCR products amplified directly from animal tissue (e.g. mouse tail) do not always resolve well when electrophoresed on agarose gels. To improve band resolution on agarose gels, add 1 µl Proteinase K to the entire 50 µl PCR sample **after** thermal cycling. You may choose not to add Proteinase K for sample types that have minimal cellular debris or protein after lysis, such as blood.
- PCR products produced by Terra Direct PCR contain 3' A-overhangs, making them compatible with T/A cloning.
- Please refer to the Terra Polymerase Mix User Manual for the Terra and General PCR troubleshooting guides.

Contact Us	
Customer Service/Ordering	Technical Support
tel: 800.662.2566 (toll-free)	tel: 800.662.2566 (toll-free)
fax: 800.424.1350 (toll-free)	fax: 800.424.1350 (toll-free)
web: takarabio.com/service	web: takarabio.com/support
e-mail: ordersUS@takarabio.com	e-mail: technical_support@takarabio.com

Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements

© 2024 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

This document has been reviewed and approved by the Quality Department.