



Map of pAcGFP1-C Vector. All restriction sites shown are unique.

### Description

pAcGFP1-C Vector is a mammalian expression vector that encodes a green fluorescent protein (GFP) from *Aequorea coerulea*. The fluorescent protein coding sequence in this construct has been human-codon-optimized for efficient expression and enhanced brightness. AcGFP1 protein has an excitation maximum at 475 nm and an emission maximum at 505 nm. Sequences flanking AcGFP1 have been converted to a Kozak consensus translation initiation site (1) to further increase translation efficiency in eukaryotic cells. A gene of interest can be added in-frame downstream of the AcGFP1 coding sequence and expressed as a fusion protein to the C-terminus of AcGFP1. SV40 polyadenylation signals downstream of the AcGFP1 gene direct proper processing of the 3' end of the mRNA transcript. In addition, the vector also contains a SV40 origin of replication in mammalian cells expressing the SV40T-antigen. A neomycin resistance cassette (Neo<sup>r</sup>) containing an SV40 early promoter, the neomycin/kanamycin resistance gene of Tn5, and polyadenylation signals from the Herpes simplex virus thymidine kinase (HSV TK) gene, allows selection of stable transformants in eukaryotic cells using G418. A bacterial promoter upstream of the gene expresses kanamycin resistance in *E. coli*. pAcGFP1-C Vector also provides a pUC origin of replication for propagation in *E. coli* and an f1 origin for single-stranded DNA production.

### Use

Provided with the Matchmaker™ Chemiluminescent CoIP Vector Set (Cat. No. 630458), the pAcGFP1-C Vector allows the generation of a N-terminal AcGFP1-bait fusion. Any bait sequence in our current Matchmaker Yeast Two Hybrid System 3 pGBKT7 bait vector can be PCR-amplified using the AcGFP BD FWD/REV universal In-Fusion primer set, and directionally and efficiently cloned into the Sall/HindIII sites of the pAcGFP1-C using our In-Fusion Dry-Down PCR Cloning Kit (Cat. No. 639602). The recombinant AcGFP1 vector can be transfected into mammalian cells using any standard transfection method. The resulting expressed AcGFP1-bait fusion protein will automatically be in-frame and can be used to monitor transfection efficiency as well as AcGFP1-bait fusion expression and localization *in vivo* since the AcGFP1-bait fusion retains the fluorescent properties of the native AcGFP1 protein.

(PR6Y2127; published 9 November 2006)



Clontech

United States/Canada  
800.662.2566

Asia Pacific  
+1.650.919.7300

Europe  
+33.(0)1.3904.6880

Japan  
+81.(0)77.543.6116

Clontech Laboratories, Inc.  
A Takara Bio Company  
1290 Terra Bella Ave.  
Mountain View, CA 94043  
Technical Support (US)  
E-mail: tech@clontech.com  
www.clontech.com

**Location of features**

- P CMV ie  
Start: 1 End: 589
- AcGFP1  
Start: 613 End: 1329
- SV40 early poly A  
Start: 1541 End: 1591
- f1 ori  
Start: 1638 End: 2093
- SV40 ori  
Start: 2434 End: 2569
- Kan R/Neo R  
Start: 2618 End: 3410
- Kanamycin/neomycin resistance gene:  
Neomycin phosphotransferase coding sequences:  
start codon (ATG): 2629–2631  
stop codon: 3421–3423  
G to A mutation to remove Pst I site: 2811  
C to A mutation to remove BssH II site: 3157
- HSVTK poly A  
Start: 3648 End: 3666
- pUC ori  
Start: 3997 End: 4640

**Propagation in *E. coli***

- Suitable host strains: DH5 $\alpha$ , HB101 and other general purpose strains. Single-stranded DNA production requires a host containing an F plasmid such as JM101 or XL1-Blue.
- Selectable marker: plasmid confers resistance to kanamycin (50  $\mu$ g/ml) in *E. coli* hosts.
- *E. coli* replication origin: pUC
- Copy number: ~500
- Plasmid incompatibility group: pMB1/ColE1

**References**

1. Kozak, M. (1987) *Nucleic Acids Res.* **15**(20):8125–8148.
2. Gorman, C. (1985) In *DNA cloning: A practical approach, Vol. II*. Ed. D.M. Glover. (IRL Press, Oxford, UK) pp. 143–190.

**Note:** The attached sequence file has been compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by Clontech. This vector has not been completely sequenced.

**Notice to Purchaser**

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, in vitro diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

**Living Colors® Products**

AcGFP1, DsRed, HcRed, AsRed, AmCyan, ZsGreen, ZsYellow and their variants:

Not-For-Profit Entities: Orders may be placed in the normal manner by contacting your local representative or Clontech Customer Service at 650.919.7300. At its discretion, Clontech grants Not-For-Profit Entities a non-exclusive, personal, limited license to use this product for non-commercial life science research use only. Such license specifically excludes the right to sell or otherwise transfer this product, its components or derivatives thereof to third parties. No modifications to the protein coding sequence may be made without express written permission from Clontech. Any other use of this product requires a license from Clontech. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at [licensing@clontech.com](mailto:licensing@clontech.com).

For-Profit Entities wishing to use this product are required to obtain a license from Clontech. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at [licensing@clontech.com](mailto:licensing@clontech.com).

This product is the subject of U.S. patents.

**CMV Sequence****ProLabel™ Products**

This product is intended to be used for research purposes only. It is not to be used for drug or diagnostic purposes nor is it intended for human use. Clontech products may not be resold, modified for resale, or used to manufacture commercial products without written approval of Clontech Laboratories, Inc.

**PCR****DISCLAIMER OF RIGHTS**

No license is conveyed with the purchase of this product under U.S. Patents Nos. 5,804,375, 5,994,056 and 6,171,785, and corresponding patent claims outside the U.S., relating to the 5' nuclease and dsDNA-binding dye processes. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc.  
Clontech is a Takara Bio Company. ©2006