#### pDsRed2-1 Vector Information

PT3602-5 Cat. No. 632405



**Restriction Map and Multiple Cloning Site (MCS) of pDsRed2-1.** Unique restriction sites are in bold. The *Not* I site is part of the DsRed2 stop codon.

#### Description

pDsRed2-1 encodes DsRed2, a DsRed variant engineered for faster maturation and lower non-specific aggregation. Derived from the Discosoma sp. red fluorescent protein (drFP583; 1), DsRed2, like its progenitor DsRed1, contains a series of silent base-pair changes that correspond to human codon-usage preferences for high expression in mammalian cells (2). In addition to these changes, DsRed2 contains six amino acid substitutions: V105A, I161T, and S197A, which result in the more rapid appearance of red fluorescence in transfected cell lines; and R2A, K5E, and K9T, which prevent the protein from aggregating. (DsRed2 may, however, form the same tetrameric structure as DsRed1 [3].) In mammalian cell cultures when DsRed2 is expressed constitutively, red-emitting cells can be detected by fluorescence microscopy within 24 hours of transfection. Large insoluble aggregates of protein, often observed in bacterial and mammalian cell systems expressing DsRed1, are dramatically reduced in cells expressing DsRed2. The faster-maturing, more soluble red fluorescent protein is also well tolerated by host cells; mammalian cell cultures transfected with DsRed2 show no obvious signs of reduced viability - in those cell lines tested, cells expressing DsRed2 display the same morphology (e.g., adherence, light-refraction) and growth characteristics as non-transfected controls.

pDsRed2-1 is a promoterless DsRed2 vector that can be used to monitor transcription from different promoters and promoter/enhancer combinations inserted into the multiple cloning site (MCS). Sequences upstream of DsRed2 have been converted to a Kozak consensus translation initiation site (4) to increase translation efficiency in eukaryotic cells. SV40 polyadenylation signals downstream of the DsRed2 gene direct proper processing of the 3' end of the DsRed2 mRNA. The vector backbone contains an SV40 origin for replication in mammalian cells expressing the SV40T antigen, a pUC origin of replication for propagation in *E. coli*, and an f1 origin for single-stranded DNA production. A neomycin-resistance cassette (Neo<sup>r</sup>) allows stably transfected eukaryotic cells to be selected using G418. This cassette consists of the SV40 early promoter, the neomycin/kanamycin resistance gene of Tn5, and polyadenylation signals from the Herpes simplex virus thymidine kinase (HSVTK) gene. A bacterial promoter upstream of the cassette expresses kanamycin resistance in *E. coli*. (PR063557; published 30 June 2010)



United States/Canada 800.662.2566 Asia Pacific +1.650.919.7300 Europe +33.(0)1.3904.6880 Japan +81.(0)77.543.6116

Clontech Laboratories, Inc. A Takara Bio Company 1290 Terra Bella Ave. Mountain View, CA 94043 Technical Support (US) E-mail: tech@clontech.com www.clontech.com

#### pDsRed2-1

# Use

DsRed2 can be used as an *in vivo* reporter of gene expression. Promoters should be cloned into the pDsRed2-1 MCS upstream from the DsRed2 coding sequence. Without addition of a functional promoter, this vector will not express DsRed2. The recombinant DsRed2 vector can be transfected into mammalian cells using any standard transfection method. If required, stable transfectants can be selected using G418 (5).

# Location of features

- MCS: 12-83
- Discosoma sp. human codon-optimized Red Fluorescent Protein (DsRed2) gene Kozak consensus translation initiation site: 90–100 Start codon (ATG): 97–99; Stop codon: 772–774
- SV40 early mRNA polyadenylation signal Polyadenylation signals: 926–931 & 955–960 mRNA 3' ends: 964 & 976
- f1 single-strand DNA origin: 1023–1478 (Packages noncoding strand of DsRed2.)
- Ampicillin resistance (b-lactamase) promoter –35 region: 1540–1545; –10 region: 1563–1568 Transcription start point: 1575
- SV40 origin of replication: 1819–1954
- SV40 early promoter Enhancer (72-bp tandem repeats): 1650–1723 & 1724–1795 21-bp repeats: 1799–1819, 1820–1840 & 1842–1862 Early promoter element: 1875–1881 Major transcription start points: 1871, 1909, 1915 & 1920
- Kanamycin/neomycin resistance gene Neomycin phosphotransferase coding sequences: Start codon (ATG): 2003–2005; stop codon: 2795–2797 G→A mutation to remove *Pst* I site: 2185 C→A (Arg→Ser) mutation to remove *Bss*H II site: 2531
- Herpes simplex virus (HSV) thymidine kinase (TK) polyadenylation signal Polyadenylation signals: 3033–3038 & 3046–3051
- pUC plasmid replication origin: 3382-4025

# **Sequencing Primer Location**

• DsRed1-C Sequencing Primer (No.632388; 5'-AGCTGGACATCACCTCCCACAACG-3'): 689-712

# Propagation in E. Coli

- Suitable host strains: DH5α, HB101 and other general purpose strains. Single-stranded DNA production requires a host containing an F plasmid such as JM109 or XL1-Blue.
- Selectable marker: plasmid confers resistance to kanamycin (50 µg/ml) to E. coli hosts.
- E. coli replication origin: pUC
- Copy number: ≈500
- Plasmid incompatibility group: pMB1/Col E1

# **Red Fluorescent Protein (DsRed2)**

• Excitation/Emission Maxima: 558 nm / 583 nm

# References

- 1. Matz, M. V., et al. (1999) Nature Biotech. 17:969–973.
- 2. Haas, J., et al. (1996) Curr. Biol. 6:315-324.
- 3. Yarbrough, D., et al. (2001) Proc. Natl. Acad. Sci. USA 98:462–467.
- 4. Kozak, M. (1987) Nucleic Acids Res. 15:8125–8148.
- 5. Gorman, C. (1985) In DNA cloning: A Practical Approach, Vol. II. Ed. D. M. Glover. (IRL Press, Oxford, U.K.), pp. 143–190.

#### **Notice to Purchaser**

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, in vitro diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

The RCFP's (including DsRedExpress and DsRedExpress2) are covered by one or more of the following U.S. Patent Nos. 7,166,444; 7,157,565; 7,217,789; 7,338,784; 7,338,783; 7,537,915 6,969,597, 7,150,979 and 7,442,522.

#### Living Colors Fluorescent Protein Products:

Not-For-Profit Entities: Orders may be placed in the normal manner by contacting your local representative or Clontech Customer Service at 650.919.7300. At its discretion, Clontech grants Not-For-Profit Entities a non-exclusive, personal, limited license to use this product for non-commercial life science research use only. Such license specifically excludes the right to sell or otherwise transfer this product, its components or derivatives thereof to third parties. No modifications to the protein coding sequence may be made without express written permission from Clontech. Any other use of this product requires a license from Clontech. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at licensing@clontech.com. For-Profit Entities wishing to use this product are required to obtain a license from Clontech. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at licensing@clontech.com.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc., unless noted otherwise. Clontech is a Takara Bio Company. ©2010 Clontech Laboratories, Inc.