

Use

pDsRed-Monomer-N1 can be used to construct fusions to the N-terminus of DsRed-Monomer. If a fusion construct retains the fluorescent properties of the native DsRed-Monomer protein, its expression can be monitored by flow cytometry and its localization *in vivo* can be determined by fluorescence microscopy. The target gene must be cloned into pDsRed-Monomer-N1 so that it is in frame with the DsRed-Monomer coding sequence, with no intervening in-frame stop codons. The inserted gene must include an initiating ATG codon. Recombinant pDsRed-Monomer-N1 can be transfected into mammalian cells using any standard transfection method. If required, stable transfectants can be selected using G418 (4). pDsRed-Monomer-N1 can also be used as a cotransfection marker; the unmodified vector will express DsRed-Monomer.

The DsRed1-N Sequencing Primer (Cat. No. 632387) can be used to sequence genes cloned adjacent to the 5' end of the DsRed-Monomer coding region.

For Western blotting, the Living Colors® DsRed Polyclonal Antibody (Cat. No. 632496) can be used to recognize the DsRed-Monomer protein. However, to generate optimal results it may be necessary to use a higher concentration of antibody than recommended on the DsRed Polyclonal Antibody Certificate of Analysis.

Location of features

- Human cytomegalovirus (CMV) immediate early promoter: 1–589
Enhancer region: 59–465; TATA box: 554–560
Transcription start point: 583
C→G mutation to remove *Sac* I site: 569
- Multiple Cloning Site: 591–671
- Human codon-optimized DsRed-Monomer gene
Kozak consensus translation initiation site: 672–682
Start codon (ATG): 679–681; Stop codon: 1354–1356
Amino acid substitutions (DsRed→DsRed-Monomer)
GCC→GAC (Ala-2 to Asp) mutation: 682–684
TCC→AAC (Ser-3 to Asn) mutation: 685–687
TCC→ACC (Ser-4 to Thr) mutation: 688–690
AAG→GAG (Lys-5 to Glu) mutation: 691–693
AAC→GAC (Asn-6 to Asp) mutation: 694–696
CGC→CAG (Arg-13 to Gln) mutation: 715–717
ACC→TCC (Thr-21 to Ser) mutation: 739–741
GAG→TAC (Glu-26 to Tyr) mutation: 754–756
CGC→AAG (Arg-36 to Lys) mutation: 784–786
CAC→ACC (His-41 to Thr) mutation: 799–801
AAC→CAG (Asn-42 to Gln) mutation: 802–804
GTG→GCC (Val-44 to Ala) mutation: 808–810
AAG→CAG (Lys-47 to Gln) mutation: 817–819
GTG→GCC (Val-71 to Ala) mutation: 889–891
AAG→ATG (Lys-83 to Met) mutation: 925–927
AAG→ACC (Lys-92 to Thr) mutation: 952–954
GTG→TCC (Val-96 to Ser) mutation: 964–966
ACC→GAG (Thr-106 to Glu) mutation: 994–996
ACC→CAG (Thr-108 to Gln) mutation: 1000–1002
TCC→ACC (Ser-117 to Thr) mutation: 1027–1029
ATC→AAG (Ile-125 to Lys) mutation: 1051–1053
TCC→GCC (Ser-131 to Ala) mutation: 1069–1071
ATG→GCC (Met-141 to Ala) mutation: 1099–1101
GCC→CCC (Ala-145 to Pro) mutation: 1111–1113
CGC→AAG (Arg-149 to Lys) mutation: 1123–1125
CGC→CAG (Arg-153 to Gln) mutation: 1135–1137
CAC→TCC (His-162 to Ser) mutation: 1162–1164
AAG→CAC (Lys-163 to His) mutation: 1165–1167
CTG→ACC (Leu-174 to Thr) mutation: 1198–1200
GTG→TGC (Val-175 to Cys) mutation: 1201–1203
GAG→GAC (Glu-176 to Asp) mutation: 1204–1206
TCC→ACC (Ser-179 to Thr) mutation: 1213–1215
ATC→GTG (Ile-180 to Val) mutation: 1216–1218
ATG→AAG (Met-182 to Lys) mutation: 1222–1224
TAC→AAC (Tyr-192 to Asn) mutation: 1252–1254

TAC→CAC (Tyr-193 to His) mutation: 1255–1257
TCC→AAC (Ser-203 to Asn) mutation: 1285–1287
ATC→GTG (Ile-210 to Val) mutation: 1306–1308
CGC→CAC (Arg-216 to His) mutation: 1324–1326
ACC→GCC (Thr-217 to Ala) mutation: 1327–1329
GGC→GCC (Gly-219 to Ala) mutation: 1333–1335
CAC→TCC (His-222 to Ser) mutation: 1342–1344
CTG→GGC (Leu-223 to Gly) mutation: 1345–1347
TTC→TCC (Phe-224 to Ser) mutation: 1348–1350
CTG→CAG (Leu-225 to Gln) mutation: 1351–1353

- SV40 early mRNA polyadenylation signal
Polyadenylation signals: 1510–1515 & 1539–1544; mRNA 3' ends: 1548 & 1560
- f1 single-strand DNA origin: 1607–2062 (Packages the noncoding strand of DsRed-Monomer)
- Bacterial promoter for expression of Kan^r gene:
–35 region: 2124–2129; –10 region: 2147–2152
Transcription start point: 2159
- SV40 origin of replication: 2403–2538
- SV40 early promoter
Enhancer (72-bp tandem repeats): 2236–2307 & 2308–2379
21-bp repeats: 2383–2403, 2404–2424 & 2426–2446
Early promoter element: 2459–2465
Major transcription start points: 2455, 2493, 2499 & 2504
- Kanamycin/neomycin resistance gene
Neomycin phosphotransferase coding sequences: Start codon (ATG): 2587–2589; Stop codon: 3379–3381
G→A mutation to remove *Pst* I site: 2769
C→A (Arg to Ser) mutation to remove *Bss*H II site: 3115
- Herpes simplex virus (HSV) thymidine kinase (TK) polyadenylation signal
Polyadenylation signals: 3617–3622 & 3630–3635
- pUC plasmid replication origin: 3966–4609

Sequencing primer location

- DsRed1-N Sequencing Primer (Cat. No. 632387; 5'-GTACTGGAAGTGGGGGGACAG-3'): 879–859

Propagation in *E. coli*

- Suitable host strains: DH5 α , HB101 and other general purpose strains. Single-stranded DNA production requires a host containing an F plasmid, such as the JM109 or XL1-Blue strains.
- Selectable marker: plasmid confers resistance to kanamycin (50 μ g/ml) in *E. coli* hosts.
- *E. coli* replication origin: pUC
- Copy number: high
- Plasmid incompatibility group: pMB1/ColE1

Excitation and emission maxima of DsRed-Monomer

- Excitation maximum = 557 nm
- Emission maximum = 592 nm

References

1. Matz, M. V., et al. (1999) *Nature Biotech.* **17**:969–973.
2. Haas, J., et al. (1996) *Curr. Biol.* **6**:315–324.
3. Kozak, M. (1987) *Nucleic Acids Res.* **15**:8125–8148.
4. Gorman, C. (1985) In *DNA cloning: A Practical Approach, Vol. II*. Ed. D. M. Glover. (IRL Press, Oxford, U.K.), pp. 143–190.

Note: The attached sequence file has been compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by Clontech Laboratories, Inc.. This vector has not been completely sequenced.

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

Not-For-Profit Entities: Orders may be placed in the normal manner by contacting your local representative or Clontech Customer Service at 650.919.7300. At its discretion, Clontech grants Not-For-Profit Entities a non-exclusive, personal, limited license to use this product for non-commercial life science research use only. Such license specifically excludes the right to sell or otherwise transfer this product, its components or derivatives thereof to third parties. No modifications to the protein coding sequence may be made without express written permission from Clontech. Any other use of this product requires a license from Clontech. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at licensing@clontech.com. For-Profit Entities wishing to use this product are required to obtain a license from Clontech. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at licensing@clontech.com or click here for more information.

The DsRed Monomer and the Fruit Fluorescent Proteins are covered by one or more of the following U.S. Patents: 7,157,566; 7,393,923; 7,005,511 and 7,250,298.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc., unless noted otherwise. Clontech is a Takara Bio Company. ©2010 Clontech Laboratories, Inc.