


pSRE-SEAP Vector Information

PT3277-5

GenBank Accession No.: Submission in progress.

Catalog No. 631901



TB=Transcription Blocker

Restriction Map of pSRE-SEAP. All restriction sites are unique.**Description:**

pSRE-SEAP is designed for monitoring the induction of the serum response element (SRE; 1) and the mitogen-activated protein (MAP) kinase signal transduction pathway (2). pSRE-SEAP contains the secreted alkaline phosphatase (SEAP) reporter gene (3–5). This vector also contains three tandem copies of the SRE consensus sequence fused to a TATA-like promoter (P_{TAL}) region from the Herpes simplex virus thymidine kinase (HSV-TK) promoter. After transcription factors bind SRE, transcription is induced and the reporter gene is activated.

The SEAP coding sequence is followed by the SV40 late polyadenylation signal to ensure proper, efficient processing of the SEAP transcript in eukaryotic cells. Located upstream of SRE is a synthetic transcription blocker (TB), which is composed of adjacent polyadenylation and transcription pause sites for reducing background transcription (6). The vector backbone also contains an f1 origin for single-stranded DNA production, a pUC origin of replication, and an ampicillin resistance gene for propagation and selection in *E. coli*.

Use:

pSRE-SEAP is designed to measure the binding of transcription factors to SRE, providing a direct measurement of activation for this pathway. For example, the addition of serum or growth factors to the cell-culture medium induces the binding of transcription factors to SRE, which initiates transcription of SEAP. Alternatively, you can cotransfect this vector with an expression vector containing your gene of interest to monitor pathway activation. The secreted SEAP enzyme can be assayed directly from the culture medium using one of Clontech's Great EscAPe Chemiluminescence Detection Kits (Cat. Nos. 631701, 631704). In addition, the SEAP assay permits time-course studies not possible with assays dependent on cell lysates. The pSRE-SEAP Vectors can be transfected into mammalian cells by any standard method. For selecting stable clones, cotransfect with a vector containing an antibiotic resistance gene, such as neomycin, hygromycin, or puromycin.

**Clontech**

United States/Canada
800.662.2566

Asia Pacific
+1.650.919.7300

Europe
+33.(0)1.3904.6880

Japan
+81.(0)77.543.6116

Clontech Laboratories, Inc.
A Takara Bio Company
1290 Terra Bella Ave.
Mountain View, CA 94043
Technical Support (US)
E-mail: tech@clontech.com
www.clontech.com

(PR93947; published 05 May 1999)

Location of features:

- Serum Response Element (SRE; 1): 27–86
- TATA-like promoter (P_{TAL}): 93–241
- Secreted alkaline phosphatase (SEAP) gene:
 - SEAP coding sequences:
 - start codon (ATG): 269–271; stop codon: 1826–1828
 - signal peptide: 269–319
 - mature protein: 320–1825
 - C-terminal extension to SEAP: 1787–1825
- SV40 late mRNA polyadenylation signal: 1939–1944
 - mRNA 3' end: 1958
- pUC plasmid replication origin: 2337–2980
- Ampicillin resistance gene:
 - Promoter: –35 region: 4058–4053; –10 region: 4035–4030
 - Transcription start point: 4023
 - Ribosome binding site: 4000–3996
 - β -lactamase coding sequences:
 - start codon (ATG): 3988–3986; stop codon: 3130–3128
 - β -lactamase signal peptide: 3988–3920
 - β -lactamase mature protein: 3919–3131
- f1 single-strand DNA origin (packages the noncoding strand of SEAP): 4120–4575
- Transcription blocker (TB): 4706–4859
 - Synthetic polyadenylation site (7): 4706–4754
 - Transcription pause site from human $\alpha 2$ globin gene (8): 4768–4859

Propagation in *E. coli*:

- Suitable host strains: DH5 α and other general purpose strains. Single-stranded DNA production requires a host containing an F' episome such as JM109.
- Selectable marker: plasmid confers resistance to ampicillin (50 μ g/ml) to *E. coli* hosts.
- *E. coli* replication origin: pUC
- Copy number: ~500
- Plasmid incompatibility group: pMB1/Col E1

References:

1. Latchman, D. S. (1995) In *Eukaryotic Transcription Factors*, 2nd Ed. (Academic Press, San Diego), pp. 72.
2. Ransone, L. & Verma, I. M. (1990) *Ann. Rev. Cell Biol.* 6:539–557.
3. Yang, T.T., et al. (July 1994) *CLONTECHniques* IX(3):1–5.
4. Berger, J., et al. (1988) *Gene* 66:1–10.
5. Cullen, B. R. & Malim, M. H. (1992) *Methods Enzymol.* 216:362–368.
6. Eggermont, J. & Proudfoot, N. (1993) *EMBO J.* 12:2539–2548.
7. Levitt, N., et al. (1989) *Genes Dev.* 3:1019–1025.
8. Enriquez-Harris, P., et al. (1991) *EMBO J.* 10:1833–1842.

Notice to Purchaser

The attached sequence file has been compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by Clontech. This vector has not been completely sequenced.

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc. Clontech is a Takara Bio Company. ©2006