

Archival nanowell sequencing

I. Introduction

Archival nanowell sequencing (Arc-well) describes high-throughput single-cell DNA-seq library preparation from nuclei isolated from fresh-frozen and FFPE tissues. The protocol was developed by scientists in Dr. Nicholas Navin's lab at the UT MD Anderson Cancer Center, and was originally published in Wang, Kaile et al. Archival single-cell genomics reveals persistent subclones during DCIS progression. *Cell* **186**, 3968–3982.e15 (2023).

[Access the paper »](#)

NOTE: This user-generated protocol is provided for general information only and is not directly supported, endorsed, or validated by Takara Bio. Takara Bio gives no warranties and makes no claims about the provided protocol. For questions, including dispense program setup and pre-experiment consultation, please contact our Field Support team at field_support@takarabio.com.

The Arc-well assay on the ICELL8® cx Single-Cell System has six dispense steps as listed below:




Figure 1. Schematic overview of Arc-well workflow on the ICELL8 cx system.

Using ICELL8 cx CELLSTUDIO® Software, set up an Arc-well application with the dispense steps shown in Figure 2. For additional information, please see the [ICELL8 cx CELLSTUDIO Application Creation Quick Start Guide](#). The software will generate a 384-well source plate layout file, which you can use as a guide when aliquoting nuclei suspension or reagent master mix into the corresponding wells.

Startup	Maintenance	SMART-Seq Pro	Arc-Well	Arc-Well Single Dispense
Chip ID: 72 x 72 : 350 nl				
Dispense Cells and Controls (35 nl)				
Scan chip				
Dispense Lysis (35 nl filtered)				
Dispense Tagmentation (35 nl filtered)				
Dispense Index 1 (35 nl filtered)				
Dispense Index 2 (35 nl filtered)				
Dispense PCR (35 nl filtered)				

	Cells	PosCtrl	NegCtrl	Lysis	Tagmentation	PCR	Index 1	Index 2																
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A	1	17	33		65	81	97	113	129	145	161		193	209	225	241	257	273						369
B	2	18	34		66	82	98	114	130	146	162		194	210	226	242	258	274						
C	3	19	35		67	83	99	115	131	147	163		195	211	227	243	259	275						
D	4	20	36		68	84	100	116	132	148	164		196	212	228	244	260	276						
E					69	85	101	117					197	213	229	245								
F					70	86	102	118					198	214	230	246								
G					71	87	103	119					199	215	231	247								
H					72	88	104	120					200	216	232	248								
I					73	89	105	121					201	217	233	249								
J					74	90	106	122					202	218	234	250								
K					75	91	107	123					203	219	235	251								
L					76	92	108	124					204	220	236	252								
M			45		77	93	109	125					205	221	237	253								
N			46		78	94	110	126					206	222	238	254								
O			47		79	95	111	127					207	223	239	255								
P			48		80	96	112	128					208	224	240	256								384

Figure 2. Arc-well dispense protocol and plate setup in the ICELL8 cx CELLSTUDIO software. Individual dispense and chip-scanning steps are shown on the left. Color-coded wells on the 384-well source plate indicate the correct position for aliquoting (right).

II. Additional Materials Required

Nuclei isolation

Reagents and equipment	Source	Catalog number
Formalin solution, neutral buffered, 10%	Millipore Sigma	HT5012-1CS
Finesse ME+	Thermo Scientific	A77500016
gentleMACS Dissociator with Heaters	Miltenyi Biotec	130-096-427
FFPE Tissue Dissociation Kit	Miltenyi Biotec	130-118-052
1.5 ml LoBind tubes	Fisher Scientific	22431021
Dimethyl sulfoxide (DMSO)	Millipore Sigma	D2650

NST-DAPI buffer for nuclei suspension

Reagents and equipment	Source	Catalog number
NaCl	Invitrogen	AM9760G
MgCl ₂	Invitrogen	AM9530G
Tris base	Fisher Scientific	BP154-1
CaCl ₂	Sigma-Aldrich	21115-1ML
BSA	Sigma-Aldrich	A2058-100G
Nonidet P-40	US Biological	N3500
DAPI	Invitrogen	D1306

Nuclei dispense

Reagents and equipment	Source	Catalog number
ICELL8 350v Chip	Takara Bio	640019
ICELL8 Blank Chip Reagent Kit (contains Second Diluent)	Takara Bio	640196
ICELL8 Loading Kit – B	Takara Bio	640206
ICELL8 Collection Kit – L	Takara Bio	640212
MSND 384-Well Source Plate and Seals (430-000025)	Takara Bio	640018
DPBS	Millipore Sigma	D8537

Lysis mix

Reagents and equipment	Source	Catalog number
Tris-HCl (pH 7.4)	Millipore Sigma	T2194
Tween-20	Andwin Scientific	NC9022994
Triton X-100	Fisher Scientific	ICN19485483
QIAGEN Protease (30 AU)	Qiagen	19157

Tagmentation mix

Reagents and equipment	Source	Catalog number
Illumina® Tagment DNA Enzyme and Buffer Large Kit	Illumina	20034198

Index 1, Index 2, and PCR mix

Reagents and equipment	Source	Catalog number
KAPA HotStart PCR Kit, with dNTPs	Roche	KK2502
0.5 M EDTA	Thermo Scientific	R1021
Nuclease-free water	Takara Bio	0912
MgCl ₂	Invitrogen	AM9530G

Please see the Appendix for additional Index information.

III. Protocol

A. Nuclei isolation from FFPE tissues

1. Prepare NST-DAPI buffer using the recipe below. Filter-sterilize the NST-DAPI buffer and store at 4°C in the dark.

Component	Volume
NST solution	800 ml
<ul style="list-style-type: none"> • 146 mM NaCl • 10 mM Tris base (pH 7.8) • 1 mM CaCl₂ • 21 mM MgCl₂ • 0.05% (w/v) BSA • 0.5% (v/v) NP-40 	
DAPI solution	200 ml
<ul style="list-style-type: none"> • 106 mM MgCl₂ • 10 mg DAPI 	
Total volume	1,000 ml

2. Section FFPE blocks to generate 1–3 50 μm scrolls using a microtome.
3. Generate single-cell suspensions from the scrolls using the gentleMACS Dissociator with Heaters and the FFPE tissue dissociation kit according to the manufacturer's recommendations with the following changes:
 - a. After the last wash with ice-cold Buffer W, transfer single-cell suspensions to 1.5 ml LoBind tubes and spin at 500g for 5 min at 4°C.
 - b. Resuspend the cell pellet in an appropriate volume of NST-DAPI buffer (200 μl–5 ml) to generate nuclei suspensions.

NOTE: Nuclei suspensions can be stored long term in 10% DMSO.
 - c. Flow-sort nuclei suspensions with the BD FACSMelody, where DAPI intensity is used to gate the desired diploid/aneuploid populations.

B. Arc-well scDNA-seq library prep on the ICELL8 cx Single-Cell System

1. Dispense nuclei and scan chip

1. Wash nuclei suspensions with 0.5X DPBS. Follow the table below and combine Second Diluent with the nuclei suspension in 0.5X DPBS for a final concentration of 28,571 nuclei per ml.

Component	Diluted stained nuclei suspension	Volume per source well (for each sample)*
Second Diluent (100X)	10 μl	1 μl
Nuclei suspension	Dilute to 28,571 nuclei per ml	Dilute to 28,571 nuclei per ml
0.5X DPBS	Up to 1,000 μl	Up to 100 μl
Total volume	1,000 μl	100 μl

*You can load up to eight different samples per chip run.

2. Using a 200 µl pipette tip, carefully load 80 µl of nuclei suspension into wells A1, A2, B1, B2, C1, C2, D1, and D2 of a 384-well source plate (blue; Figure 2). If applicable, prepare positive and negative control samples, and refer to Figure 2 to aliquot to corresponding source wells.
3. Load the 350v chip in the chip nest, place the 384-well source plate in the ICELL8 cx Single-Cell System with the A1 corner positioned at the top-right corner of the plate nest.
4. Click [Dispense cells and Controls (35 nl)].
5. Make sure the chip and source plate are properly inserted, remove the seals, and click [Done] to confirm.
6. When the dispense is completed, remove the chip for blotting, and seal the chip with RC Film using a film-sealing roller.
7. Centrifuge the chip at 1,000g for 5 min at 4°C.
8. Load the chip, remove the RC Film, and click [Scan chip].
9. Follow steps in the [ICELL8 cx Single-Cell System User Manual](#) (Section X.E, Page 62) to acquire imaging data (make sure to scan the chip under the DAPI channel), perform automatic data analysis, and save the result in the ICELL8 cx CellSelect® software.

2. Dispense Lysis mix

1. Prepare Lysis mix as below:

Component	Volume
Lysis buffer	180 µl
<ul style="list-style-type: none"> • 30 mM Tris-HCl, pH 8.0 • 5% Tween • 0.5% Triton X-100 	
Protease (1.36 AU/ml)	20 µl
Total volume	200 µl

4. Aliquot 50 µl of freshly prepared Lysis mix into wells A3, B3, C3, and D3 (orange; Figure 2) of a 384-well source plate and seal the plate.
5. Load chip and source plate in the correct orientation. Remove imaging film from the chip and sealing film from the source plate.
6. Click [Dispense Lysis (35 nl filtered)].
7. When the dispense is completed, remove the chip for blotting, and seal the chip with RC Film using a film-sealing roller.
8. Centrifuge the chip at 1,000g for 5 min at 4°C.
9. Put the chip onto the ICELL8 cx thermal cycler and run the Lysis program below:

59.7°C	5 sec
54.5°C	30 min
79°C	11 sec
75.3°C	15 min
4°C	hold

3. Dispense Tagmentation mix

1. Prepare Tagmentation mix as below:

Component	Volume
2x TD buffer	144 μ l
TDE1, Illumina	16 μ l
Total volume	160 μ l

2. Aliquot 40 μ l of freshly prepared Tagmentation mix into wells A11, B11, C11, and D11 (yellow; Figure 2) of a 384-well source plate and seal the plate.
3. Load chip and source plate in the correct orientation. Remove RC Film from the chip and sealing film from the source plate.
4. Click [Dispense Tagmentation (35 nl filtered)].
5. When the dispense is completed, remove the chip for blotting, and seal the chip with RC Film using a film-sealing roller.
6. Centrifuge the chip at 3,220g for 5 min at 4°C.
7. Put the chip onto the ICELL8 cx thermal cycler and run the Tagmentation program below:

59.7°C	5 sec
54.5°C	<ul style="list-style-type: none"> • 8 min (for FFPE sample) • 12 min (for fresh sample)

4. Dispense Index 1 (with Neutralization mix)

1. Prepare Neutralization master mix as below:

Component	Volume in master mix (for 72 Index wells)
5x Kapa Hifi Fidelity buffer	850 μ l
25 mM dNTP	136 μ l
0.5 M EDTA	85 μ l
Nuclease-free water	544 μ l
Total volume	1,615 μ l

2. Aliquot 19 μ l of the Neutralization master mix into each well for wells A5–P8, A9–D9, and A10–D10 (coral; Figure 2) on the source plate. Spike in 1 μ l of Index 1 and seal:

Component	Per well in wells A5–P8, A9–D9, and A10–D10
Neutralization master mix	19 μ l
100 μ M Index 1	1 μ l
Total volume per Index well	20 μ l

3. Centrifuge the source plate and the chip (with RC Film on) at 3,220g (minimum 2,600g) for 3 min at 4°C.
4. Load chip and source plate in the correct orientation. Remove RC Film from the chip and sealing film from the source plate.
5. Click [Dispense Index 1 (35 nl filtered)].
6. When the dispense is completed, remove the chip for blotting, and seal the chip with RC Film using a film-sealing roller.
7. Centrifuge the chip at 3,220g for 5 min at 4°C.
8. Put the chip onto the ICELL8 cx thermal cycler and run the Neutralization program below:

54.9°C	5 sec
49.4°C	30 min
4°C	hold

5. Dispense Index 2

1. Prepare master mix as below:

Component	Volume in master mix (for 72 Index wells)
5x Kapa Hifi Fidelity buffer	850 µl
1 M MgCl ₂	29.58 µl
Nuclease-free water	735.42 µl
Total volume	1,615 µl

2. Aliquot 19 µl of the above master mix into each well for wells A13–P16, A17–D17, and A18–D18 (cyan; Figure 2) on the source plate. Spike in 1 µl of Index 2 and seal:

Component	Per well in wells A13–P16, A17–D17, and A18–D18
Master mix	19 µl
100 µM Index 2	1 µl
Total volume per Index well	20 µl

3. Centrifuge the source plate and the chip (with RC Film on) at 3,220g (minimum 2,600g) for 3 min at 4°C.
4. Load chip and source plate in the correct orientation. Remove RC Film from the chip and sealing film from the source plate.
5. Click [Dispense Index 2 (35 nl filtered)].
6. When the dispense is completed, remove the chip for blotting, and seal the chip with RC Film using a film-sealing roller.
7. Centrifuge the chip at 3,220g for 5 min at 4°C.

6. Dispense PCR mix

1. Prepare PCR mix as below:

Component	Volume
5x Kapa HiFi Fidelity buffer	40 μ l
KAPA HiFi HotStart DNA Polymerase (1 U/ μ l)	40 μ l
Nuclease-free water	120 μ l
Total volume	200 μ l

2. Aliquot 50 μ l of freshly prepared PCR mix into wells M3, N3, O3, and P3 (pink; Figure 2) of a 384-well source plate and seal the plate.
3. Load chip and source plate in the correct orientation. Remove RC film from the chip and sealing film from the source plate.
4. Click [Dispense PCR (35 nl filtered)].
5. When the dispense is completed, remove the chip for blotting, and seal the chip with RC Film using a film-sealing roller.
6. Centrifuge the chip at 3,220g for 5 min at 4°C.
7. Put the chip onto the ICELL8 cx thermal cycler and run the PCR program below:

72.1°C 8 min

99.6°C 30 sec

X cycles*:

99.6°C 20 sec

57.5°C 5 sec

62.7°C 30 sec

72.1°C 1 min

72.1°C 2 min

4°C hold

*10-12 cycles for frozen samples; 14-16 cycles for FFPE samples

8. Centrifuge the chip at 3,220g for 5 min at 4°C.

7. Extract, purify, and QC the library

1. Open ICELL8 Collection Kit – L and label the Collection Tube with the engraved chip number.
2. Assemble the collection module by attaching the Collection Tube to the Collection Fixture.
3. Carefully peel the film from the chip.
4. With the nanowells facing down, place the chip into the assembled collection module.
5. Seal the chip and the top of the collection module with a Collection Film.
6. Using a balance or blank chip, assemble another collection module.
7. Centrifuge both collection modules at 3,220g (minimum 2,600g) for 10 min at 4°C.

8. Perform 1.8X bead purification with AmPure XP beads on 50 µl of the library.
9. QC the library with Agilent BioAnalyzer 2100. Use the bioanalyzer results to determine library quality and average size for qPCR.
10. Quantitate the library using the Library Quantification Kit (Cat. No. 638324). Refer to the [Library Quantification Kit User Manual](#) for instructions. Use the average size as determined by the bioanalyzer to calculate the molar library concentration.
11. Use the qPCR results from the Library Quantification Kit to determine the library quantity for sequencing.
12. The library can be sequenced on Illumina NextSeq® 2000 or NovaSeq™ 6000 with 8 bp dual indexing sequencing at a target of 1 million reads per cell.

IV. Appendix

Index 1		
Well position	Name	Sequence
A5	ArcS501	AATGATACGGCGACCACCGAGATCTACACTAGATCGCTCGTCGGCAGCGTC
B5	ArcS502	AATGATACGGCGACCACCGAGATCTACACTATCCTCTTCGTCGGCAGCGTC
C5	ArcS503	AATGATACGGCGACCACCGAGATCTACACAGAGTAGATCGTCGGCAGCGTC
D5	ArcS504	AATGATACGGCGACCACCGAGATCTACACGTAAGGAGTCGTCGGCAGCGTC
E5	ArcS505	AATGATACGGCGACCACCGAGATCTACACACTGCATATCGTCGGCAGCGTC
F5	ArcS506	AATGATACGGCGACCACCGAGATCTACACAAGGAGTATCGTCGGCAGCGTC
G5	ArcS507	AATGATACGGCGACCACCGAGATCTACACCTAAGCCTTCGTCGGCAGCGTC
H5	ArcS508	AATGATACGGCGACCACCGAGATCTACACCGTCTAATTCGTCGGCAGCGTC
I5	ArcS509	AATGATACGGCGACCACCGAGATCTACACTCTCTCCGTCGTCGGCAGCGTC
J5	ArcS510	AATGATACGGCGACCACCGAGATCTACACTCGACTAGTCGTCGGCAGCGTC
K5	ArcS511	AATGATACGGCGACCACCGAGATCTACACTTCTAGCTTCGTCGGCAGCGTC
L5	ArcS512	AATGATACGGCGACCACCGAGATCTACACCCTAGAGTTCGTCGGCAGCGTC
M5	ArcS513	AATGATACGGCGACCACCGAGATCTACACGCGTAAGATCGTCGGCAGCGTC
N5	ArcS514	AATGATACGGCGACCACCGAGATCTACACCTATTAAGTCGTCGGCAGCGTC
O5	ArcS515	AATGATACGGCGACCACCGAGATCTACACAAGGCTATTCGTCGGCAGCGTC
P5	ArcS516	AATGATACGGCGACCACCGAGATCTACACGAGCCTTATCGTCGGCAGCGTC
A6	ArcS517	AATGATACGGCGACCACCGAGATCTACACTTATGCGATCGTCGGCAGCGTC
B6	ArcS518	AATGATACGGCGACCACCGAGATCTACACTAAGCGTTTCGTCGGCAGCGTC
C6	ArcS519	AATGATACGGCGACCACCGAGATCTACACTCCGTCTTTCGTCGGCAGCGTC
D6	ArcS520	AATGATACGGCGACCACCGAGATCTACACTTCTGTGTTTCGTCGGCAGCGTC
E6	ArcS521	AATGATACGGCGACCACCGAGATCTACACTCTGCTGTTTCGTCGGCAGCGTC

F6	ArcS522	AATGATACGGCGACCACCGAGATCTACA CTGGAGGTTTCGTCGGCAGCGTC
G6	ArcS523	AATGATACGGCGACCACCGAGATCTACA CTGATACGTTTCGTCGGCAGCGTC
H6	ArcS524	AATGATACGGCGACCACCGAGATCTACA CTGCATAGTTTCGTCGGCAGCGTC
I6	ArcS525	AATGATACGGCGACCACCGAGATCTACA CTGCGATCTTCGTCGGCAGCGTC
J6	ArcS526	AATGATACGGCGACCACCGAGATCTACA CTTCCTGCTTCGTCGGCAGCGTC
K6	ArcS527	AATGATACGGCGACCACCGAGATCTACA CTAGTGACTTCGTCGGCAGCGTC
L6	ArcS528	AATGATACGGCGACCACCGAGATCTACA CTACAGGATTCGTCGGCAGCGTC
M6	ArcS529	AATGATACGGCGACCACCGAGATCTACA CTGTGGTTGTCGTCGGCAGCGTC
N6	ArcS530	AATGATACGGCGACCACCGAGATCTACA CTACTAGTCTCGTCGGCAGCGTC
O6	ArcS531	AATGATACGGCGACCACCGAGATCTACA CTCGAAGTGTCGTCGGCAGCGTC
P6	ArcS532	AATGATACGGCGACCACCGAGATCTACA CTAAACGCTGTCGTCGGCAGCGTC
A7	ArcS533	AATGATACGGCGACCACCGAGATCTACA CTTGGTATGTCGTCGGCAGCGTC
B7	ArcS534	AATGATACGGCGACCACCGAGATCTACA CTGAACTGGTTCGTCGGCAGCGTC
C7	ArcS535	AATGATACGGCGACCACCGAGATCTACA CTACTTTCGGTTCGTCGGCAGCGTC
D7	ArcS536	AATGATACGGCGACCACCGAGATCTACA CTCTCACGGTTCGTCGGCAGCGTC
E7	ArcS537	AATGATACGGCGACCACCGAGATCTACA CTACAGACGTGTCGTCGGCAGCGTC
F7	ArcS538	AATGATACGGCGACCACCGAGATCTACA CTTGCCTAATCGTCGGCAGCGTC
G7	ArcS539	AATGATACGGCGACCACCGAGATCTACA CTTTAACATCGTCGGCAGCGTC
H7	ArcS540	AATGATACGGCGACCACCGAGATCTACA CTACCGTAGACCTCGTCGGCAGCGTC
I7	ArcS541	AATGATACGGCGACCACCGAGATCTACA CTATTTGCGTTCGTCGGCAGCGTC
J7	ArcS542	AATGATACGGCGACCACCGAGATCTACA CTACATCCAGGATCGTCGGCAGCGTC
K7	ArcS543	AATGATACGGCGACCACCGAGATCTACA CTCTGGCGATCGTCGGCAGCGTC
L7	ArcS544	AATGATACGGCGACCACCGAGATCTACA CAATCTACATCGTCGGCAGCGTC
M7	ArcS545	AATGATACGGCGACCACCGAGATCTACA CCGATAGGGTTCGTCGGCAGCGTC
N7	ArcS546	AATGATACGGCGACCACCGAGATCTACA CCGTGAAGGTCGTCGGCAGCGTC
O7	ArcS547	AATGATACGGCGACCACCGAGATCTACA CATCGAATGTCGTCGGCAGCGTC
P7	ArcS548	AATGATACGGCGACCACCGAGATCTACA CTCAAGAGCTCGTCGGCAGCGTC
A8	ArcS549	AATGATACGGCGACCACCGAGATCTACA CCGCCACGTTTCGTCGGCAGCGTC
B8	ArcS550	AATGATACGGCGACCACCGAGATCTACA CCCCCTTGGATCGTCGGCAGCGTC
C8	ArcS551	AATGATACGGCGACCACCGAGATCTACA CATTACCGTTTCGTCGGCAGCGTC
D8	ArcS552	AATGATACGGCGACCACCGAGATCTACA CAGTCCGAGTCGTCGGCAGCGTC
E8	ArcS553	AATGATACGGCGACCACCGAGATCTACA CACTTGTTGTCGTCGGCAGCGTC

F8	ArcS554	AATGATACGGCGACCACCGAGATCTACACGTAATACATCGTCGGCAGCGTC
G8	ArcS555	AATGATACGGCGACCACCGAGATCTACACGGCGTCTATCGTCGGCAGCGTC
H8	ArcS556	AATGATACGGCGACCACCGAGATCTACACGCGCTGCTTCGTCGGCAGCGTC
I8	ArcS557	AATGATACGGCGACCACCGAGATCTACACGTGCCATTTTCGTCGGCAGCGTC
J8	ArcS558	AATGATACGGCGACCACCGAGATCTACACAACACCTATCGTCGGCAGCGTC
K8	ArcS559	AATGATACGGCGACCACCGAGATCTACACCTCCGAACTCGTCGGCAGCGTC
L8	ArcS560	AATGATACGGCGACCACCGAGATCTACACCAACGGCATCGTCGGCAGCGTC
M8	ArcS561	AATGATACGGCGACCACCGAGATCTACACCAATGTAGTCGTCGGCAGCGTC
N8	ArcS562	AATGATACGGCGACCACCGAGATCTACACGGCTACCCTCGTCGGCAGCGTC
O8	ArcS563	AATGATACGGCGACCACCGAGATCTACACAAAGTCCGTCGTCGGCAGCGTC
P8	ArcS564	AATGATACGGCGACCACCGAGATCTACACTTCCGCGGTCGTCGGCAGCGTC
A9	ArcS565	AATGATACGGCGACCACCGAGATCTACACAGGCACTTTTCGTCGGCAGCGTC
B9	ArcS566	AATGATACGGCGACCACCGAGATCTACACCTTCAGTGTCGTCGGCAGCGTC
C9	ArcS567	AATGATACGGCGACCACCGAGATCTACACGCCGGTAGTCGTCGGCAGCGTC
D9	ArcS568	AATGATACGGCGACCACCGAGATCTACACTTCAATCCTCGTCGGCAGCGTC
A10	ArcS569	AATGATACGGCGACCACCGAGATCTACACCCACACACTCGTCGGCAGCGTC
B10	ArcS570	AATGATACGGCGACCACCGAGATCTACACATATTATCTCGTCGGCAGCGTC
C10	ArcS571	AATGATACGGCGACCACCGAGATCTACACCCGAAGCATCGTCGGCAGCGTC
D10	ArcS572	AATGATACGGCGACCACCGAGATCTACACGTATCGGTTTCGTCGGCAGCGTC

Index 2		
Well position	Name	Sequence
A13	ArcN701	CAAGCAGAAGACGGCATAACGAGATTGCCTTAGTCTCGTGGGCTCGG
B13	ArcN702	CAAGCAGAAGACGGCATAACGAGATCTAGTACGGTCTCGTGGGCTCGG
C13	ArcN703	CAAGCAGAAGACGGCATAACGAGATTTCTGCCTGTCTCGTGGGCTCGG
D13	ArcN704	CAAGCAGAAGACGGCATAACGAGATGCTCAGGAGTCTCGTGGGCTCGG
E13	ArcN705	CAAGCAGAAGACGGCATAACGAGATAGGAGTCCGTCTCGTGGGCTCGG
F13	ArcN706	CAAGCAGAAGACGGCATAACGAGATCATGCCTAGTCTCGTGGGCTCGG
G13	ArcN707	CAAGCAGAAGACGGCATAACGAGATGTAGAGAGGTCTCGTGGGCTCGG
H13	ArcN708	CAAGCAGAAGACGGCATAACGAGATCCTCTCTGGTCTCGTGGGCTCGG
I13	ArcN709	CAAGCAGAAGACGGCATAACGAGATAGCGTAGCGTCTCGTGGGCTCGG
J13	ArcN710	CAAGCAGAAGACGGCATAACGAGATCAGCCTCGGTCTCGTGGGCTCGG

K13	ArcN711	CAAGCAGAAGACGGCATAACGAGATTGCCTCTTGTCTCGTGGGCTCGG
L13	ArcN712	CAAGCAGAAGACGGCATAACGAGATTCTCTACGTCTCGTGGGCTCGG
M13	ArcN713	CAAGCAGAAGACGGCATAACGAGATTCATGAGCGTCTCGTGGGCTCGG
N13	ArcN714	CAAGCAGAAGACGGCATAACGAGATCCTGAGATGTCTCGTGGGCTCGG
O13	ArcN715	CAAGCAGAAGACGGCATAACGAGATTAGCGAGTGTCTCGTGGGCTCGG
P13	ArcN716	CAAGCAGAAGACGGCATAACGAGATGTAGCTCCGTCTCGTGGGCTCGG
A14	ArcN717	CAAGCAGAAGACGGCATAACGAGATTACTACGCGTCTCGTGGGCTCGG
B14	ArcN718	CAAGCAGAAGACGGCATAACGAGATAGGCTCCGGTCTCGTGGGCTCGG
C14	ArcN719	CAAGCAGAAGACGGCATAACGAGATGCAGCGTAGTCTCGTGGGCTCGG
D14	ArcN720	CAAGCAGAAGACGGCATAACGAGATCTGCGCATGTCTCGTGGGCTCGG
E14	ArcN721	CAAGCAGAAGACGGCATAACGAGATGAGCGCTAGTCTCGTGGGCTCGG
F14	ArcN722	CAAGCAGAAGACGGCATAACGAGATACTGATCGGTCTCGTGGGCTCGG
G14	ArcN723	CAAGCAGAAGACGGCATAACGAGATTAGCTGCAGTCTCGTGGGCTCGG
H14	ArcN724	CAAGCAGAAGACGGCATAACGAGATGACGTCGAGTCTCGTGGGCTCGG
I14	ArcN725	CAAGCAGAAGACGGCATAACGAGATATCACGTTGTCTCGTGGGCTCGG
J14	ArcN726	CAAGCAGAAGACGGCATAACGAGATCGATGTTTGTCTCGTGGGCTCGG
K14	ArcN727	CAAGCAGAAGACGGCATAACGAGATTTAGGCATGTCTCGTGGGCTCGG
L14	ArcN728	CAAGCAGAAGACGGCATAACGAGATTGACCACTGTCTCGTGGGCTCGG
M14	ArcN729	CAAGCAGAAGACGGCATAACGAGATACAGTGGTGTCTCGTGGGCTCGG
N14	ArcN730	CAAGCAGAAGACGGCATAACGAGATGCCAATGTGTCTCGTGGGCTCGG
O14	ArcN731	CAAGCAGAAGACGGCATAACGAGATCAGATCTGGTCTCGTGGGCTCGG
P14	ArcN732	CAAGCAGAAGACGGCATAACGAGATACTTGATGGTCTCGTGGGCTCGG
A15	ArcN733	CAAGCAGAAGACGGCATAACGAGATGATCAGCGGTCTCGTGGGCTCGG
B15	ArcN734	CAAGCAGAAGACGGCATAACGAGATGGCTACAGGTCTCGTGGGCTCGG
C15	ArcN735	CAAGCAGAAGACGGCATAACGAGATTGGTTGTTGTCTCGTGGGCTCGG
D15	ArcN736	CAAGCAGAAGACGGCATAACGAGATTCTCGGTTGTCTCGTGGGCTCGG
E15	ArcN737	CAAGCAGAAGACGGCATAACGAGATTAATAAGAGTCTCGTGGGCTCGG
F15	ArcN738	CAAGCAGAAGACGGCATAACGAGATACTAAGTCGTCTCGTGGGCTCGG
G15	ArcN739	CAAGCAGAAGACGGCATAACGAGATGCTGGTCTGTCTCGTGGGCTCGG
H15	ArcN740	CAAGCAGAAGACGGCATAACGAGATCTGTATTTGTCTCGTGGGCTCGG
I15	ArcN741	CAAGCAGAAGACGGCATAACGAGATTTTCATAAGTCTCGTGGGCTCGG
J15	ArcN742	CAAGCAGAAGACGGCATAACGAGATGACCCAAGGTCTCGTGGGCTCGG

K15	ArcN743	CAAGCAGAAGACGGCATAACGAGATTTATTTGGGTCTCGTGGGCTCGG
L15	ArcN744	CAAGCAGAAGACGGCATAACGAGATTTTAACGCGTCTCGTGGGCTCGG
M15	ArcN745	CAAGCAGAAGACGGCATAACGAGATCTTACTCCGTCTCGTGGGCTCGG
N15	ArcN746	CAAGCAGAAGACGGCATAACGAGATGGGAACCGGTCTCGTGGGCTCGG
O15	ArcN747	CAAGCAGAAGACGGCATAACGAGATGCATTAAGTCTCGTGGGCTCGG
P15	ArcN748	CAAGCAGAAGACGGCATAACGAGATGACCGTTTGTCTCGTGGGCTCGG
A16	ArcN749	CAAGCAGAAGACGGCATAACGAGATTTTGGATCGTCTCGTGGGCTCGG
B16	ArcN750	CAAGCAGAAGACGGCATAACGAGATATCATCATGTCTCGTGGGCTCGG
C16	ArcN751	CAAGCAGAAGACGGCATAACGAGATCGTGTTGGGTCTCGTGGGCTCGG
D16	ArcN752	CAAGCAGAAGACGGCATAACGAGATTGTTGTTAGTCTCGTGGGCTCGG
E16	ArcN753	CAAGCAGAAGACGGCATAACGAGATGGTTTACCGTCTCGTGGGCTCGG
F16	ArcN754	CAAGCAGAAGACGGCATAACGAGATGGTCGATGGTCTCGTGGGCTCGG
G16	ArcN755	CAAGCAGAAGACGGCATAACGAGATGTTCCCATGTCTCGTGGGCTCGG
H16	ArcN756	CAAGCAGAAGACGGCATAACGAGATATTGGCCGGTCTCGTGGGCTCGG
I16	ArcN757	CAAGCAGAAGACGGCATAACGAGATTCATTCCCGTCTCGTGGGCTCGG
J16	ArcN758	CAAGCAGAAGACGGCATAACGAGATCCGAATACGTCTCGTGGGCTCGG
K16	ArcN759	CAAGCAGAAGACGGCATAACGAGATATAGCTGAGTCTCGTGGGCTCGG
L16	ArcN760	CAAGCAGAAGACGGCATAACGAGATAGATAAATGTCTCGTGGGCTCGG
M16	ArcN761	CAAGCAGAAGACGGCATAACGAGATATCTCGGGGTCTCGTGGGCTCGG
N16	ArcN762	CAAGCAGAAGACGGCATAACGAGATAGACATTAGTCTCGTGGGCTCGG
O16	ArcN763	CAAGCAGAAGACGGCATAACGAGATGAATTGGCGTCTCGTGGGCTCGG
P16	ArcN764	CAAGCAGAAGACGGCATAACGAGATGCACGGCGGTCTCGTGGGCTCGG
A17	ArcN765	CAAGCAGAAGACGGCATAACGAGATTCGGTCAGGTCTCGTGGGCTCGG
B17	ArcN766	CAAGCAGAAGACGGCATAACGAGATTCGAAATGGTCTCGTGGGCTCGG
C17	ArcN767	CAAGCAGAAGACGGCATAACGAGATTGGCAAGCGTCTCGTGGGCTCGG
D17	ArcN768	CAAGCAGAAGACGGCATAACGAGATTGGTAGAAGTCTCGTGGGCTCGG
A18	ArcN769	CAAGCAGAAGACGGCATAACGAGATCGTCACGTGTCTCGTGGGCTCGG
B18	ArcN770	CAAGCAGAAGACGGCATAACGAGATCGCGACAGTCTCGTGGGCTCGG
C18	ArcN771	CAAGCAGAAGACGGCATAACGAGATAAGTTTAAGTCTCGTGGGCTCGG
D18	ArcN772	CAAGCAGAAGACGGCATAACGAGATGTTGTGGTGTCTCGTGGGCTCGG

Contact Us	
Customer Service/Ordering	Technical Support
tel: 800.662.2566 (toll-free)	tel: 800.662.2566 (toll-free)
fax: 800.424.1350 (toll-free)	fax: 800.424.1350 (toll-free)
web: takarabio.com/service	web: takarabio.com/support
e-mail: ordersUS@takarabio.com	e-mail: technical_support@takarabio.com

Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product’s web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

© 2024 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

United States/Canada: +1.800.662.2566 • Asia Pacific: +1.650.919.7300 • Europe: +33.(0)1.3904.6880 • Japan: +81.(0)77.565.6999

1.24 US