

Scaling up plate-based single-cell profiling: introducing automated STORM-seq



Cedric Uyttingco¹, Chris Waddling³, Kamila Koprowska³, Arnab Mandal², Smriti Chawla², Mike Covington¹, Matthieu Pesant¹, Benjamin K. Johnson⁴, Mary F. Rhodes⁴, Saloni Pasta¹, Anitha Jayaprakash¹, Michael Phelan¹, Yue Yun¹, Hui Shen⁴ & Andrew Farmer¹

1. Takara Bio USA, Inc., 2560 Orchard Pkwy, San Jose, CA, 95131, USA
 2. DSS Takara Bio India, Pvt. Ltd., Jayanagar 4th Block, Bengaluru, KA, 560011, India
 3. SPT Labtech, Melbourn Science Park, Melbourn Hertfordshire SG8 6HB, United Kingdom
 4. Department of Epigenetics, Van Andel Institute, Grand Rapids, MI, 49503, USA

*Corresponding Authors

Background

Full-length, single-cell RNA-seq (scRNA-seq) allows for more complete transcriptomic profiling compared to 5' or 3' end-counting methods. However, labor-intensive, non-modular workflows with high reagent costs limit scalability. Advancements in laboratory automation, protocol design, and reaction miniaturization have overcome some of these constraints. In collaboration with the Van Andel Institute, we present Single-cell Total RNA Miniaturized sequencing (STORM-seq), which pairs Takara Bio's SMART[™] technology with SPT Labtech's mosquito[®] HV genomics liquid handling platform. This miniaturized, plate-based scRNA-seq approach is scalable, automated, and offers a solution for single-cell transcriptomic profiling (i.e., coding/noncoding transcripts, isoforms, gene fusions, etc.) at reduced cost and increased cell throughput.

Methods

Single viable K562 cells, peripheral blood mononuclear cells (PBMCs), or renal cell carcinoma cells (RCCs) were sorted into 384-well plates using fluorescence-activated cell sorting (FACS). Following lysis and RNA fragmentation, RNA was reverse transcribed into cDNA. Unique molecular identifiers (UMIs) and unique dual indices (UDIs) were incorporated into the cDNA library following a modified protocol from Johnson and Rhodes et al. (2022, bioRxiv). Miniaturized (1/6) reactions were facilitated using mosquito HV genomics (SPT Labtech). After indexing, pooling, and cleanup, the cDNA libraries were treated with the ZapR[™] ribo-depletion technology to remove rRNA-derived cDNAs. The libraries were amplified and sequenced using an Illumina[®] NextSeq[®] 500. Sequencing results were analyzed using Cogent[™] NGS Software (Takara Bio).

1 STORM-seq with UMIs: a flexible and modular workflow




Figure 1. STORM-seq with UMIs features an easy workflow that generates Illumina-compatible scRNA-seq libraries. Actual processing time is around 8 hours, and it may vary depending on the number of samples, third-party instrumentation, and cycling conditions. The workflow includes multiple safe stopping points: post-cell sorting, after pooling (Step 4), and after PCR2 amplification (Step 7).

Takara Bio USA, Inc.
 United States/Canada: +1.800.662.2566 • Asia Pacific: +1.850.919.7300 • Europe: +33.(0)1.3904.6880 • Japan: +81.(0)77.565.6999
 FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES. © 2024 Takara Bio Inc. All Rights Reserved.
 All trademarks are the property of Takara Bio Inc. or its affiliates in the U.S. and/or other countries or their respective owners.
 Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com

2 Comparison of single-cell library recovery using a manual or automated approach




Figure 2. 96-well manual and 384-well automated STORM-seq reaction plates exhibit >95% recovery of libraries and single cells. Total barcode reads (top) and % mapped reads (bottom) were determined for both the 96-well manual and the 384-well automated STORM-seq protocol. Total barcode reads, and % mapped reads were largely uniform across the reaction plate for both the manual and automated protocol. Negative controls are shown in wells G12 and H12 for the 96-well manual protocol and A1 and P24 for the 384-well automated protocol.

3 Sequencing metrics, UMI and gene detection, and gene body coverage

Sequencing metrics		
RNA source	Single K562 cells	
Plate type	384-well	
Final concentration, total yield	137 ng/ul, 2.189 ng	
Replicates	N = 78	
Mean	CV (%)	
Mapped reads (%)	98	2
Uniquely mapped reads (%)	90	2
Exonic reads (%)	76	11
Intronic reads (%)	19	41.4
Intergenic reads (%)	5	25
Number of UMIs	53,444	9
Number of genes	10,248	12
Number of transcripts	17,126	16




Figure 3. Performance of libraries generated using the automated STORM-seq workflow. cDNA libraries were generated from single K562 cells and sequenced on an Illumina NextSeq instrument (2 x 75 bp). Sequences were normalized to 1 x 10⁶ paired-end reads. Panel A. Sequencing metrics. The libraries also showed >98% strand specificity (data not shown). Panel B. UMI and gene saturation plots. TPM=transcripts per million. Panel C. Gene body coverage plot.

4 Gene expression is reproducible between single K562 cells




Figure 4. Gene expression data is reproducible between cells in the automated STORM-seq protocol. Panel A. Correlation plots for four randomly selected K562 cell replicates demonstrate >0.9 Pearson correlation between single-cell samples. Panel B. Heatmap showing Pearson's Correlation of gene expression derived from 381 cells.

5 Detection of multiple RNA biotypes




Figure 5. STORM-seq identifies multiple RNA biotypes. Counts of selected biotypes for the automated STORM-seq protocol identified from 381 K562 cells.

6 Detection of gene fusions in K562 cells




Figure 6. STORM-seq identifies gene fusion events in K562 cells. Panel A. Circos plot visualization and quantification of junction fusion events. Panel B. Circos plot visualization and quantification of spanning fusion events.

7 Detection of distinct single-cell populations




Figure 7. STORM-seq identifies distinct cell populations. Panel A. Heatmap with dendrogram illustrating differentially expressed genes in PBMCs. Panel B. uMAP plot from RCCs, K562 cells, and human brain control RNA. Panel C. Heatmap with dendrogram illustrating differentially expressed genes in RCCs, K562 cells, and human brain control RNA.

Conclusions

- STORM-seq is a flexible (automated or manual) and cost-efficient method for transcriptome-wide gene expression analysis of single cells with relatively even gene body coverage.
- Scalable from several to thousands of cells, STORM-seq has an adaptable experimental design to accommodate various inputs, such as cryopreserved or FACS-isolated cells from cell lines, primary cells, and tumor cells.
- Primer design and uniform gene body coverage in STORM-seq enable the detection of coding and noncoding RNA biotypes as well as different transcript isoforms and gene fusions.

References

Johnson, B. and Rhodes, M. et al. Single-cell Total RNA Miniaturized sequencing (STORM-seq) reveals differentiation trajectories of primary human fallopian tube epithelium. *bioRxiv* (2022).

Figure 1 was created with BioRender.com



For further inquiries, please email
 Saloni.Pasta@takarabio.com

800.662.2566
 takarabio.com