

# Lenti-X<sup>™</sup> Tet-Off<sup>®</sup> Advanced Vector Set

Catalog No.	Amount	Lot Number
632158 (Not sold separately)	10 µg	Specified on product label.
Sold as a part of 632163		

# Description

The Lenti-X Tet-Off Advanced Vector Set contains lentiviral expression vectors that are designed to establish an inducible expression system in mammalian cells. Lentivirus produced from the pLVX-Tet-Off Advanced Vector expresses a tetracycline-dependent transcriptional activator (tTA-Advanced) which specifically induces transcription of your gene of interest from the  $P_{\text{Tight}}$  inducible promoter in pLVX-Tight-Puro. In cells cotransduced by both lentiviral vectors, your cDNA is expressed to high levels in the absence of doxycycline or after its withdrawal. Using these vectors with the Lenti-X HTX Packaging System (Cat. No. 631247) will produce very high-titer supernatants from 293T packaging cells. The VSV-G pseudotyped lentiviruses so generated can be used to transduce virtually any mammalian cell type.

## **Package Contents**

- 10 µg pLVX-Tet-Off Advanced Vector (500 ng/µl)
- 10 µg pLVX-Tight-Puro Vector (500 ng/µl)
- 10 µg pLVX-Tight-Puro-Luc Vector (500 ng/µl)

## Storage Conditions

• Store at –20°C.

## Shelf Life

• 1 year from date of receipt under proper storage conditions.

## Storage Buffer

- 10 mM Tris-HCl (pH 8.0)
- 1 mM EDTA (pH 8.0)

## Concentration

500 ng/μl

## **Shipping Conditions**

• Dry ice  $(-70^{\circ}C)$ 

## **Product Documents**

Documents for our products are available for download at <u>takarabio.com/manuals</u> The following documents apply to this product:

- pLVX-Tet-Off Advanced Vector Information Packet
- pLVX-Tight-Puro Vector Information Packet
- pLVX-Tight-Puro-Luc Vector Information Packet

## Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA U.S. Technical Support: <u>techUS@takarabio.com</u>

# Certificate of Analysis

Lenti-X Tet-Off Advanced Vector Set (Not sold separately)

# **Quality Control Data**

# **Plasmid Identity & Purity**

• The identity of each plasmid was verified by sequencing critical elements, and by electrophoresis on a 0.8% agarose/EtBr gel after restriction enzyme digestion.:

Plasmid Name	Enzyme	Fragments (bp)
pLVX-Tet-Off Advanced	XhoI KpnI	11,578 9,718 & 1,860
pLVX-Tight-Puro	XhoI SalI	7,191 5,606 & 2,185
pLVX-Tight-Puro-Luc	XhoI SalI	9,424 7,239 & 2,185

• The purity of each plasmid was also checked by determining the  $A_{260}/A_{280}$ .

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.



# Lenti-X<sup>TM</sup> Tet-Off® Advanced Vector Set

## CATALOG NO.

632158

## NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at <u>http://www.takarabio.com</u>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

## **STATEMENT 42**

Use of the Tetracycline controllable expression systems (the "Tet Technology") is covered by a series of patents including U.S. Patent # 7541446, # 8383364, # 9181556, European patents EP # 1200607, # 1954811, #2352833 and corresponding patent claims outside these regions which are proprietary to TET Systems GmbH & Co. KG. Academic research institutions are granted an automatic license with the purchase of this product to use the Tet Technology only for internal, academic research purposes, which license specifically excludes the right to sell, or otherwise transfer, the Tet Technology or its component parts to third parties. Notwithstanding the above, academic and not-for profit research institutions whose research using the Tet Technology is sponsored by for profit organizations, which shall receive ownership to any data and results stemming from the sponsored research, shall need a commercial license agreement from TET Systems in order to use the Tet Technology. In accepting this license, all users acknowledge that the Tet Technology is experimental in nature. TET Systems GmbH & Co. KG makes no warranties, express or implied or of any kind, and hereby disclaims any warranties, representations, or guarantees of any kind as to the Tet Technology, patents, or products. All others are invited to request a license from TET Systems GmbH & Co. KG prior to purchasing these reagents or using them for any purpose. Takara Bio USA, Inc. is required by its licensing agreement to submit a report of all purchasers of the Tet-controllable expression system to TET Systems.

For license information, please contact: GSF/CEO TET Systems GmbH & Co. KG, Im Neuenheimer Feld 582 69120 Heidelberg Germany Tel: +49 6221 5880400 Fax: +49 6221 5880404 email: info@tetsystems.com

Takara Bio USA, Inc.				
1290 Terra Bella Avenue, M	Mountain View, CA 940	)43, USA		
U.S. Technical Support: tec	hUS@takarabio.com			4/17/2020
United States/Canada	Asia Pacific	Europe	Japan	4/1//2020
800.662.2566	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.565.6999	



or use the electronic licensing request form via http://www.tetsystems.com/ip-licensing/licensing/for-profit-research

### **STATEMENT 55**

cPPT Element. This product and its use are the subject to one or more of the following U.S. Pat. Nos. 8,093,042 and 10,407,695 and foreign equivalents. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot disclose information, sell or otherwise transfer this product or its components or materials made using this product or its components for any commercial purposes. If the buyer is not willing to accept the limitations of this limited use statement, Takara Bio USA, Inc. is willing to accept return of the product with a full refund. For information on purchasing a license to the DNA-Flap technology for purposes other than research, contact the Transfer of Technology Office, Institut Pasteur, 28 rue du Docteur Roux, 75 724 Paris Cedex 15 (www.pasteur.fr).

### **TRADEMARKS:**

### ©2018 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

Takara Bio USA, Inc.				
1290 Terra Bella Avenue, M	Mountain View, CA 940	)43, USA		
U.S. Technical Support: tec	hUS@takarabio.com			111512020
United States/Canada	Asia Pacific	Europe	Japan	4/17/2020
800.662.2566	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.565.6999	