

pAdenoX-ZsGreen1 (Linear)

Catalog No(s). 632261 (Not sold separately). Sold as part of 632267.	Amount 10 µl pAdenoX-ZsGreen1 (Linear) Vector (200 ng/µl)	Lot Number Specified on product label.
---	---	--

pAdenoX-ZsGreen1 (Linear) Vector Information




Figure 1. pAdenoX-ZsGreen1 (Linear) vector map.




Figure 2. pAdenoX-ZsGreen1 (Linear) In-Fusion cloning site.

Location of Features

- SV40 polyA signals: 106–903
- Screening Primer (reverse) [complementary]: 936–955
- $P_{CMV\ IE}$ (human cytomegalovirus immediate early promoter): 25353–25941
- ZsGreen1: 25966–26658
- SV40 polyA signals: 26814–26864
- ITR (inverted terminal repeat): 31825–31884
- Amp^r (ampicillin resistance gene; β -lactamase): 32785–33645
- pUC origin of replication: 34290–34463
- ITR (inverted terminal repeat): 34874–34933
- Screening Primer (forward): 34977–35001
- Ψ (packaging signal): 35066–35214
- $P_{CMV\ IE}$ (human cytomegalovirus immediate early promoter): 35254–35843

NOTICE TO PURCHASER:

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, in vitro diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with the licensing requirements listed below and described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.