
Vector Map

Catalog No.



pGFPuv Vector



U.S. Technical Support: <u>tech@clontech.com</u>

United States/Canada	Asia Pacific	Europe	Japan
800.662.2566	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.543.6116
(121913)			

Vector Map

pGFPuv Vector

Location of Features

- P_{lac} (*lac* promoter): 95–178
- *lacZ*-GFPuv fusion protein expressed in *E. coli*: 217–1005
- 5' MCS: 234–281
- GFPuv gene: 289–1005
 GFP chromophore: 481–489
 wt GFP cDNA sequence fragment (1): 289–454
 Synthetic GFP gene with "cycle 3" mutations from pBAD-GFPuv (2): 455-1007
- 3' MCS: 1007–1091
- Amp^R promoter: 1432-1536
- Amp^R (ampicillin resistance gene): 1537–2397
- ori (high-copy-number ColE1/pMB1/pBR322/pUC origin of replication): 2545–3188

Recommended Primer Location

• 331–352

References

- 1. Prasher, D. C., et al. (1992) Gene 111:229-233.
- 2. Crameri, A., et al. (1996) Nature Biotechnol. 14:315–319.

Notice to Purchaser

Clontech® products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with the licensing requirements described on the product's web page at <u>http://www.clontech.com</u>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

Clontech and the Clontech logo are trademarks of Clontech Laboratories, Inc. All other marks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech is a Takara Bio Company. ©2014 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.