


pZFHD1-2

Catalog No(s).

635082 (Not sold separately).
Sold as part of 635081.

pZFHD1-2 Vector Information



	AgeI	Apal	XmaI	BamHI	SpeI	NotI
462	ACC GGT	GGG CCC	CCC GGG	GGA TCC	ACT AGT	TGA TCA
	TGG CCA	CCC GGG	GGG CCC	CCT AGG	TGA TCA	ACT AGT
						CGC CGG CGC

pZFHD1-2 vector map and multiple cloning sites.

Clontech Laboratories, Inc.

A Takara Bio Company
1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: tech@clontech.com

United States/Canada 800.662.2566 Last update: 071612	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.543.6116
---	---------------------------------	------------------------------	-----------------------------

Location of Features

- P_{zI-1} (zI-1 promoter): 88–438
12 ZFHD1 binding sites: 88–297
 P_{IL2} min (IL2 minimal promoter): 319–438
- MCS (Multiple cloning site): 462–505
- SV40 polyA signal (SV40 polyadenylation signal): 650-655 and 679-684
- SV40 early promoter and enhancer: 688-694
- f1 ori (f1 origin of replication): 747–1202
- SV40 ori (SV40 origin of replication): 1543-1681
- Kan^r/Neo^r (Kanamycin/Neomycin resistance gene): 1727–2521
- HSV TK polyA signal (herpes simplex virus thymidine kinase polyA signal): 2757-2762 and 2770-2775
- pUC ori (pUC-derived origin of replication): 3106–3749
- Transcription blocker: 3779–3932

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with the licensing requirements described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

Clontech and the Clontech logo are trademarks of Clontech Laboratories, Inc. All other marks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech is a Takara Bio Company. ©2012 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.