

Product Components List

SMART-Seq® mRNA HT

Catalog No.
634791

Amount
4 x 96 rxns

Description

SMART-Seq mRNA HT is designed for synthesis of high-quality cDNA directly from 1–100 intact cells or ultra-low amounts of total RNA (10 pg–1 ng). By integrating locked nucleic acid (LNA) and SMART® (Switching Mechanism at 5' end of RNA Template) technologies, this kit offers remarkable sensitivity, enriches for full-length transcripts, and maintains true representation of the original mRNA pool, which are critical factors for transcriptome analysis. The simplified workflow of this kit allows users to process a large number of samples with ease. Sequencing libraries from cDNA generated with this kit have been tested for compatibility with Illumina® sequencing platforms.

Package Contents

- 4 x SMART-Seq mRNA HT (Cat. No. 634796) >> [View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for our products are available for download at takarabio.com/manuals
The following documents apply to this product:

- SMART-Seq mRNA HT User Manual

Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

© 2023 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

This document has been reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 945131, USA

U.S. Technical Support: technical_support@takarabio.com

United States/Canada
800.662.2566
(050423)

Asia Pacific
+1.650.919.7300

Europe
+33.(0)1.3904.6880

Japan
+81.(0)77.565.6999