

Product Components List

SmartChip ND™ Real-Time PCR System

Catalog No.
640290

Amount
Each

Description

The SmartChip ND Real-Time PCR System offers high-throughput genotyping and gene expression analysis using SmartChip® MyDesign chips that contain 5,184 individual nanowells. Once the SmartChip ND dispenser has dispensed the samples and reagents, the SmartChip MyDesign chip can be directly transferred to the SmartChip ND Real-Time PCR Cycler for PCR amplification.

Package Contents

- SmartChip ND Crate Components (Cat. No. 640292) (Not sold separately) >> [View Components](#)
- SmartChip ND Real-Time PCR Cycler (Cat. No. 640023) (Not sold separately) >> [View Components](#)
- SmartChip ND Real-Time PCR System Installation Kit (Cat. No. 640294) (Not sold separately) >> [View Components](#)
- SmartChip Real-Time PCR System Training Kit - Genotyping (Cat. No. 640154) (Not sold separately) >> [View Components](#)
- SmartChip Real-Time PCR System Training Kit - Gene Expression (Cat. No. 640155) (Not sold separately) >> [View Components](#)

For storage conditions, please see the Certificate of Analysis supplied with each component.

Product Documents

Documents for our products are available for download at takarabio.com/manuals
The following documents apply to this product:

- SmartChip ND and qPCRStudio™ Software User Manual
- SmartChip ND Real-Time PCR Cycler and Software User Manual

Notice to Purchaser

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product's web page at takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

© 2024 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions. Additional product, intellectual property, and restricted use information is available at takarabio.com.

This document has been reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA
U.S. Technical Support: technical_support@takarabio.com

United States/Canada	Asia Pacific	Europe	Japan
800.662.2566 (092724)	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.565.6999