

Single-Cell 250f Chip

Catalog No.

640193

Also sold as part of Cat. No. 640282 & 640188

Amount

1 chip

Lot Number

Specified on product label.

Description

The Single-Cell 250f Chip is a 250 nl volume, flat-bottom chip which enables users to test cell dispense and/or cell imaging using the Shasta™ Single Cell System (Cat. No. 640282) or ICELL8® cx Single-Cell System (Cat. No. 640188). All wells in the 5,184 nanowell chip are unprinted and have a maximum working volume of 200 nl and a depth of 1.600 mm.

Package Contents

- 1 Single-Cell 250f Chip

Storage Conditions

- Store at room temperature.

Shelf Life

- 1 year from date of receipt under proper storage conditions.

Shipping Conditions

- Room temperature

Product Documents

Documents for our products are available for download at takarabio.com/manuals

The following documents apply to this product:

- Shasta Single Cell System User Manual
- Shasta CellSelect Software User Manual
- ICELL8 cx Single-Cell System User Manual
- ICELL8 cx CellSelect Software User Manual

Quality Control Data

Each chip is tested for critical dimensions and contaminants.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: technical_support@takarabio.comUnited States/Canada
800.662.2566
(012925)Asia Pacific
+1.650.919.7300Europe
+33.(0)1.3904.6880Japan
+81.(0)77.565.6999

Single-Cell 250f Chip

CATALOG NO.

640193

NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at <http://www.takarabio.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

TRADEMARKS:

©2025 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: technical_support@takarabio.com

United States/Canada

800.662.2566

Asia Pacific

+1.650.919.7300

Europe

+33.(0)1.3904.6880

Japan

+81.(0)77.565.6999

1/29/2025