

Adeno-X™ GoStix™

Catalog No. 632270	Amount 20 tests	Lot Number Specified on product label.
------------------------------	---------------------------	--

Description

Adeno-X GoStix are designed to instantly confirm the presence of adenovirus in cell culture supernatants or lysates. The test detects adenoviral hexon proteins and can be used to determine if infectious virus was produced post-transfection, or for selecting the best time to harvest your virus. The included Hexon Control provides confirmation of GoStix test function.

Package Contents

- 20 each Adeno-X GoStix (1 test each)
- 2 vials Adeno-X Chase Buffer (3 ml/vial)
- 1 vial Hexon Control (1 test/vial)

Storage Conditions

- Store all components at room temperature upon receipt.
- Reconstitute the Hexon Control in 100 µl of Adeno-X Chase Buffer. Use 100 µl per cassette. Use immediately—DO NOT store.

Shelf Life

- 1 year from date of receipt under proper storage conditions.

Shipping Conditions

- Room temperature

Product Documents

Documents for our products are available for download at takarabio.com/manuals

The following documents apply to this product:

- Adeno-X GoStix Protocol-At-A-Glance

Quality Control Data

Adeno-X GoStix cassettes were tested with 0 and 5 ng of purified adenoviral hexon proteins. The intensity of the test band (T) was as expected for the hexon-containing sample, while a distinct control band (C) was observed in both tests. When GoStix cassettes were tested with the resuspended Hexon Control according to the protocol, distinct test and control bands were visible within 5–20 min of adding Adeno-X Chase Buffer.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: technical_support@takarabio.com

United States/Canada 800.662.2566 (080724)	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.565.6999
--	---------------------------------	------------------------------	-----------------------------

Adeno-X™ GoStix™

CATALOG NO.

632270

NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at <http://www.takarabio.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

STATEMENT 89

This Product is protected by one or more patents from the family consisting of: US9945850 and any corresponding patents, divisionals, continuations, patent applications and foreign filings sharing common priority with the same family.

TRADEMARKS:

©2024 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: technical_support@takarabio.com

United States/Canada

800.662.2566

Asia Pacific

+1.650.919.7300

Europe

+33.(0)1.3904.6880

Japan

+81.(0)77.565.6999

8/7/2024