Certificate of Analysis



Adeno-X™ Promoterless Vector Set (with DsRed-Express)

Catalog No. Amount Lot Number

632259 (not sold separately; sold as a 10 rxns Specified on product label. part of 632265)

Description

The Adeno-X Promoterless Vector Set (with DsRed-Express) is supplied with the Adeno-X Adenoviral System 3 (Universal, Red) [Cat. No. 632265], which lets you place an expression cassette of your choice in an adenoviral vector format. The system provides a perfect platform for the tissue-specific expression of your gene of interest, or for the efficient delivery of shRNA or microRNA. pAdenoX-PRLS-DsRedExpress is a prelinearized, adenoviral vector that is ready for the insertion of an expression cassette via In-Fusion® HD PCR cloning technology. Simply PCR-amplify the expression cassette and combine it with pAdenoX-PRLS-DsRedExpress in an In-Fusion HD Cloning reaction. In-Fusion HD Cloning is fast, simple, precise, and efficient, making Adeno-X Adenoviral System 3 the most advanced, commercially-available, adenoviral gene delivery tool. In addition, independent, constitutive expression of the fluorescent protein DsRed-Express from the E3 region of the adenoviral genome allows fluorescence visualization of infected cells during rescue, amplification, and target cell infection.

Package Contents

- 10 μl pAdenoX-PRLS-DsRedExpress (Linear) Vector (200 ng/μl)
- 50 μl Adeno-X Screening Primer Mix 3 (10 μM)
- 20 μl Adeno-X Control Fragment (50 ng/μl)

Storage Conditions

- Store at -20° C.
- Spin briefly to recover contents.
- Avoid repeated freeze/thaw cycles.

Shelf Life

1 year from date of receipt under proper storage conditions.

Shipping Conditions

• Dry ice (-70°C)

Product Documents

Documents for our products are available for download at <u>takarabio.com/manuals</u> The following documents apply to this product:

- Adeno-X Adenoviral System 3 User Manual (PT5177-1)
- pAdenoX-PRLS-DsRedExpress (Linear) Vector Map

Quality Control Data

Plasmid Identity & Purity

Digestion with the indicated restriction enzymes produced fragments of the indicated sizes on a 0.8% agarose/EtBr gel:

Vector	Enzymes	Fragment Sizes
pAdenoX-PRLS-DsRedExpress	NheI	687, 1508, 2827, 3848, 10,342, & 15,147 bp
	XhoI	595, 1445, 2466, 5653, 9700 & 14,500 bp
	PacI	10, 2965, & 31,380 bp

- Vector identity was confirmed by sequencing.
- A_{260}/A_{280} : 1.8–2.0

Functional Testing

The Adeno-X Adenoviral System 3 (Universal, red) was tested using the control fragment (*lacZ*) according to the protocol described in the Adeno-X Adenoviral System 3 User Manual (PT5177-1). Chemically competent Stellar *E. coli* cells were transformed with 1.5 μl of the In-Fusion reaction mixture. After 60 min at 37°C in SOC medium, the cells were plated on agar containing 100 μg/ml ampicillin. Transformants were grown at 37°C for 24–30 hrs. PCR colony screening with the Adeno-X Screening Primer Mix revealed that >50% of the resultant colonies contained recombinant adenoviral DNA.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

(013120) Page 2 of 2



Adeno-XTM Promoterless Vector Set (with DsRed-Express)

CATALOG NO.

632259

NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at http://www.takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

STATEMENT 25

This product is covered by U.S. Patent No. 6,303,362.

STATEMENT 72

Living Colors Fluorescent Protein Products:

Not-For-Profit Entities: Orders may be placed in the normal manner by contacting your local representative or Takara Bio USA, Inc. Customer Service. Any and all uses of this product will be subject to the terms and conditions of the Non-Commercial Use License Agreement (the "Non-Commercial License"), a copy of which can be found below. As a condition of sale of this product to you, and prior to using this product, you must agree to the terms and conditions of the Non-Commercial License. Under the Non-Commercial License, Takara Bio USA, Inc. grants Not-For-Profit Entities a non-exclusive, non-transferable, non-sublicensable and limited license to use this product for internal, noncommercial scientific research use only. Such license specifically excludes the right to sell or otherwise transfer this product, its components or derivatives thereof to third parties. No modifications to the product may be made without express written permission from Takara Bio USA, Inc. Any other use of this product requires a different license from Takara Bio USA, Inc. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at licensing@takarabio.com.

For-Profit Entities wishing to use this product are required to obtain a license from Takara Bio USA, Inc. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at licensing@takarabio.com.

Not-For-Profit Non-Commercial Use License:

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA U.S. Technical Support: techUS@takarabio.com

United States/Canada **Asia Pacific** Europe Japan 800.662.2566

4/6/2021

Notice to Purchaser



STATEMENT 24

DsRed2, DsRedExpress, DsRedExpress2, E2-Crimson, HcRed and pTimer are covered by one or more of the following U.S. Patent Nos. 7,166,444; 7,150,979; 7,157,565; 7,183,399; 7,858,844; 8,093,450; 8,461,769 and non-U.S. international equivalents in some countries.

TRADEMARKS:

©2011 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA U.S. Technical Support: techUS@takarabio.com

 United States/Canada
 Asia Pacific
 Europe
 Japan

 800.662.2566
 +1.650.919.7300
 +33.(0)1.3904.6880
 +81.(0)77.565.6999