

Certificate of Analysis

Illumina® Indexing Primer Set

Catalog Nos.

634844 (Not sold separately; sold as a part of 634836 & 634837)

634845 (Not sold separately; sold as a part of 634838 & 634839)

Amount

12 rxns

48 rxns

Lot Number

Specified on product label.

Specified on product label.

Description

The Illumina Indexing Primer Set contains PCR primers for the amplification of indexed, paired-end Illumina-compatible sequencing libraries.

Package Contents

634844 (12 rxns)	634845 (48 rxns)	
12 µl	48 µl	Universal Forward PCR Primer (12.5 µM)
12 µl	48 µl	Reverse PCR Primer Index 1 (12.5 µM)
12 µl	48 µl	Reverse PCR Primer Index 2 (12.5 µM)
12 µl	48 µl	Reverse PCR Primer Index 3 (12.5 µM)
12 µl	48 µl	Reverse PCR Primer Index 4 (12.5 µM)
12 µl	48 µl	Reverse PCR Primer Index 5 (12.5 µM)
12 µl	48 µl	Reverse PCR Primer Index 6 (12.5 µM)
12 µl	48 µl	Reverse PCR Primer Index 7 (12.5 µM)
12 µl	48 µl	Reverse PCR Primer Index 8 (12.5 µM)
12 µl	48 µl	Reverse PCR Primer Index 9 (12.5 µM)
12 µl	48 µl	Reverse PCR Primer Index 10 (12.5 µM)
12 µl	48 µl	Reverse PCR Primer Index 11 (12.5 µM)
12 µl	48 µl	Reverse PCR Primer Index 12 (12.5 µM)

Storage Conditions

- Store at –20°C

Expiration Date

- Specified on product label.

Shipping Conditions

- Dry ice

Product Documents

Documents for our products are available for download at takarabio.com/manuals

The following documents apply to this product:

- SMARTer® Stranded RNA-Seq Kit User Manual

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: technical_support@takarabio.com

United States/Canada
800.662.2566
(011623)

Asia Pacific
+1.650.919.7300

Europe
+33.(0)1.3904.6880

Japan
+81.(0)77.565.6999

Certificate of Analysis

Cat. Nos. 634844 & 634845

Illumina Indexing Primer Set (Not sold separately)

Sold as part of 634836, 634837, 634838 & 634839

Quality Control Data

All primers were verified by ESI analysis to contain a single major peak with a measured molecular weight within 10 Da of the predicted molecular weight.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

Illumina® Indexing Primer Set

CATALOG NOS.

634844 & 634845

NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at <http://www.takarabio.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

TRADEMARKS:

©2023 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: technical_support@takarabio.com

United States/Canada

800.662.2566

Asia Pacific

+1.650.919.7300

Europe

+33.(0)1.3904.6880

Japan

+81.(0)77.565.6999

1/16/2023