Certificate of Analysis



pLVX-DD-AcGFP1-Actin Vector

Catalog No. **Lot Number Amount** $10 \mu g$

631077 (Not sold separately) Sold as a part of 631076.

Specified on product label.

Description

This lentiviral expression vector encodes the ProteoTunerTM destabilization domain (DD), fused to the genes for human cytoplasmic β-actin and a green fluorescent protein derived from Aequorea coerulescens (AcGFP1). Adding or removing the DD's stabilization ligand precisely regulates the amount of stabilized DD-AcGFP1-Actin that is available to be incorporated into actin filaments. Thus, adding and removing Shield1 allows for a "pulse-chase" like experimental design, which enables you to monitor actin filament network dynamics.

To package the vector into high-titer, replication-incompetent lentiviral particles, we recommend using the Lenti-XTM HT Packaging System (Cat. No. 632160) and the Lenti-X 293T Cell Line (Cat. No. 632180). The resulting lentivirus can then be used to transduce virtually any mammalian cell type.

Package Contents

10 μg pLVX-DD-AcGFP1-Actin Vector

Storage Conditions

- Store at -20° C.
- Spin briefly to recover contents.
- Avoid repeated freeze/thaw cycles.

Shelf Life

1 year from date of receipt under proper storage conditions.

Storage Buffer

10 mM Tris-HCl (pH 8.0), 1 mM EDTA (pH 8.0)

Concentration

 $500 \text{ ng/}\mu\text{l}$

Shipping Conditions

Dry ice $(-70^{\circ}C)$

Product Documents

Documents for our products are available for download at takarabio.com/manuals The following documents apply to this product:

- Lenti-X Actin Dynamics Monitoring Kit Protocol-At-A-Glance
- pLVX-DD-AcGFP1-Actin Vector Information

Certificate of Analysis

pLVX-DD-AcGFP1-Actin Vector (Not sold separately)

Sold as a part of Cat. No. 631076

Cat. No. 631077

Propagation in *E. coli*

- Recommended host strain: StellarTM Competent Cells (Cat. No. 636763).
- Selectable marker: Plasmids confer resistance to ampicillin (100 μg/ml) in *E. coli* hosts.
- E. coli replication origin: pUC

Quality Control Data

Plasmid Identity & Purity

Digestion with the indicated restriction enzymes produced fragments of the indicated sizes on a 0.8% agarose/EtBr gel:

Vector	Enzyme(s)	Fragment(s)
pLVX-DD-AcGFP1-Actin Vector	EcoRI	10.22 kb
	Acc65I	6 47 & 3 75 kb

- Vector identity was confirmed by sequencing.
- A₂₆₀/A₂₈₀: 1.8–2.0

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

(112718) Page 2 of 2



pLVX-DD-AcGFP1-Actin Vector

CATALOG NO.

631077

NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at http://www.takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

STATEMENT 39

AcGFP is covered by U.S. Patent No. 7,432,053.

STATEMENT 57

This product is covered by U.S. Patent No. 8,173,792.

STATEMENT 72

Living Colors Fluorescent Protein Products:

Not-For-Profit Entities: Orders may be placed in the normal manner by contacting your local representative or Takara Bio USA, Inc. Customer Service. Any and all uses of this product will be subject to the terms and conditions of the Non-Commercial Use License Agreement (the "Non-Commercial License"), a copy of which can be found below. As a condition of sale of this product to you, and prior to using this product, you must agree to the terms and conditions of the Non-Commercial License. Under the Non-Commercial License, Takara Bio USA, Inc. grants Not-For-Profit Entities a non-exclusive, non-transferable, non-sublicensable and limited license to use this product for internal, non-commercial scientific research use only. Such license specifically excludes the right to sell or otherwise transfer this product, its components or derivatives thereof to third parties. No modifications to the product may be made without express written permission from Takara Bio USA, Inc. Any other use of this product requires a different license from Takara Bio USA, Inc. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at licensing@takarabio.com.

For-Profit Entities wishing to use this product are required to obtain a license from Takara Bio USA, Inc. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at licensing@takarabio.com.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA U.S. Technical Support: techUS@takarabio.com

 United States/Canada
 Asia Pacific
 Europe
 Japan

 800.662.2566
 +1.650.919.7300
 +33.(0)1.3904.6880
 +81.(0)77.565.6999

11/15/2019

Notice to Purchaser



Not-For-Profit Non-Commercial Use License:

STATEMENT 55

cPPT Element. This product and its use are the subject to one or more of the following U.S. Pat. 8.093,042 and foreign equivalents. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot disclose information, sell or otherwise transfer this product, its components or materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for any commercial purposes. If the buyer is not willing to accept the limitations of this limited use statement, Takara Bio USA, Inc. is willing to accept return of the product with a full refund. For information on purchasing a license to the DNA-Flap technology for purposes other than research, contact the Transfer of Technology Office, Institut Pasteur, 28 rue du Docteur Roux, 75 724 Paris Cedex 15 (www.pasteur.fr).

TRADEMARKS:

© 2018 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA U.S. Technical Support: techUS@takarabio.com

United States/Canada Asia Pacific **Europe** Japan +1.650.919.7300 +33.(0)1.3904.6880