# Certificate of Analysis



### **Primer Mix Kit**

Catalog No. Amount Lot Number

638326 (Not sold separately) 500 rxns Specified on product label. Sold as a part of 638324

## **Description**

The Primer Mix Kit includes components needed for the Library Quantification Kit, a kit designed for the qPCR-based quantification of DNA libraries prior to Illumina® next-generation sequencing. The components in this kit include primers that recognize the Illumina adapter sequences, EASY Dilution Buffer for diluting library samples, and RNase-free water for preparing reaction mixes.

### **Package Contents**

- 2 x 1 ml 5X Primer Mix
- 3 x 1 ml EASY Dilution Buffer
- 3 x 1 ml RNase-Free Water

# **Storage Conditions**

• Store at –20°C.

### **Expiration Date**

• Specified on product label.

### Shipping Conditions

Dry ice

#### **Product Documents**

Documents for our products are available for download at <u>takarabio.com/manuals</u> The following documents apply to this product:

• Library Quantification Kit User Manual

# **Quality Control Data**

Real-time qPCR was performed using the Terra<sup>TM</sup> qPCR Direct TB Green® Premix (Takara Bio, Cat. No. 638319), DNA Standards for Library Quantification (Takara Bio, Cat. No. 638325), and the components of this kit. The resulting standard curve had an R<sup>2</sup>>0.980, and the amplification efficiency was between 0.9 and 1.1.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.



4/24/2023

# **Primer Mix Kit**

### CATALOG NO.

638326

#### NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at <a href="http://www.takarabio.com">http://www.takarabio.com</a>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

#### **TRADEMARKS:**

#### ©2023 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA U.S. Technical Support: technical\_support@takarabio.com

United States/Canada Asia Pacific Europe Japan

800.662.2566 +1.650.919.7300 +33.(0)1.3904.6880 +81.(0)77.565.6999