

Retro-X™ Tet-On® 3G Vector Set

Catalog No.	Amount	Lot Number
631192 (Not sold separately). Sold as a part of 631188.	Each	Specified on product label.

Description

When used with the Retro-X Tet-On 3G Inducible Expression System, the Retro-X Tet-On 3G vector set allows for the retroviral delivery and inducible expression of your gene of interest in a variety of mammalian cells. Target cells that express the Tet-On 3G transactivator protein and contain a gene of interest under the control of a pTRE3GV promoter will express high levels of your gene, only when cultured in the presence of doxycycline.

Package Contents

- 20 µl pRetroX-TRE3G Vector (500 ng/µl)
- 20 µl pRetroX-TRE3G-Luc Control Vector (500 ng/µl)
- 20 µl pRetroX-Tet3G Vector (500 ng/µl)

Storage Conditions

- Store plasmids at -20°C.
- Spin briefly to recover contents.
- Avoid repeated freeze/thaw cycles.

Shelf Life

- 1 year from date of receipt under proper storage conditions.

Storage Buffer

- 10 mM Tris-HCl (pH 8.0), 1 mM EDTA (pH 8.0)

Concentration

- 500 ng/µl

Shipping Conditions

- Dry ice (-70°C)

Product Documents

Documents for our products are available for download at takarabio.com/manuals

The following documents apply to this product:

- Retro-X Tet-On 3G Inducible Expression System User Manual
- pRetroX-TRE3G Vector Information
- pRetroX-TRE3G-Luc Control Vector Information
- pRetroX-Tet3G Vector Information

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: techUS@takarabio.com

United States/Canada 800.662.2566 (050418)	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.565.6999
--	---------------------------------	------------------------------	-----------------------------

Certificate of Analysis

Cat. No. 631192

Retro-X Tet-On 3G Vector Set (Not sold separately; sold as a part of 631188)

Propagation in *E. coli*

- Recommended host strain: Stellar™ Competent Cells (Cat. No. 636763).
- Selectable marker: Plasmids confer resistance to ampicillin (100 µg/ml) in *E. coli* hosts.
- *E. coli* replication origin: pUC

Quality Control Data

Plasmid Identity & Purity

- Digestion with the indicated restriction enzymes produced fragments of the indicated sizes on a 0.8% agarose/EtBr gel:

Vector	Enzyme(s)	Fragment(s)
pRetroX-TRE3G-Luc Control	BgIII	8.2 kb
pRetroX-TRE3G	XhoI	5.0 & 3.2 kb
	EcoRV	6.6 kb
	HindIII	4.8 & 1.8 kb
pRetroX-Tet3G	BamHI	8.1 kb
	SpeI	7.3 & 0.9 kb

- Vector identity was confirmed by sequencing.
- A_{260}/A_{280} : 1.8–2.0

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

Retro-X™ Tet-On® 3G Vector Set

CATALOG NO.

631192

NOTICE TO PURCHASER:

Our products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, in vitro diagnostic purposes, therapeutics, or in humans. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Takara Bio USA, Inc.

Your use of this product is also subject to compliance with the licensing requirements listed below and described on the product's web page at <http://www.takarabio.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

STATEMENT 8

For Research Use Only, Not for Use in Humans.

STATEMENT 42

Use of the Tetracycline controllable expression systems (the "Tet Technology") is covered by a series of patents including U.S. Patent # 7541446, # 8383364, # 9181556, European patents EP # 1200607, # 1954811, #2352833 and corresponding patent claims outside these regions which are proprietary to TET Systems GmbH & Co. KG. Academic research institutions are granted an automatic license with the purchase of this product to use the Tet Technology only for internal, academic research purposes, which license specifically excludes the right to sell, or otherwise transfer, the Tet Technology or its component parts to third parties. Notwithstanding the above, academic and not-for profit research institutions whose research using the Tet Technology is sponsored by for profit organizations, which shall receive ownership to any data and results stemming from the sponsored research, shall need a commercial license agreement from TET Systems in order to use the Tet Technology. In accepting this license, all users acknowledge that the Tet Technology is experimental in nature. TET Systems GmbH & Co. KG makes no warranties, express or implied or of any kind, and hereby disclaims any warranties, representations, or guarantees of any kind as to the Tet Technology, patents, or products. All others are invited to request a license from TET Systems GmbH & Co. KG prior to purchasing these reagents or using them for any purpose. Takara Bio USA, Inc. is required by its licensing agreement to submit a report of all purchasers of the Tet-controllable expression system to TET Systems.

For license information, please contact:

GSF/CEO

TET Systems GmbH & Co. KG,

Im Neuenheimer Feld 582

69120 Heidelberg

Germany

Tel: +49 6221 5880400

Fax: +49 6221 5880404

email: info@tetsystems.com

or use the electronic licensing request form via <http://www.tetsystems.com/ip-licensing/licensing/for-profit-research>

TRADEMARKS:

©2015 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: techUS@takarabio.com

United States/Canada

Asia Pacific

Europe

Japan

5/4/2018

800.662.2566

+1.650.919.7300

+33.(0)1.3904.6880

+81.(0)77.565.6999