

Yeastmaker™ Carrier DNA

Catalog No.
630440

Amount
Each

Lot Number
Specified on product label.

Description

Yeastmaker Carrier DNA is a sheared, denatured salmon sperm DNA ready for use as yeast transformation carrier DNA.

Package Contents

- 5 x 1 ml Yeastmaker Carrier DNA (10 mg/ml)

Storage Conditions

- Store at -20°C.

Expiration Date

- Specified on product label.

Storage Buffer

- 1X TE

Shipping Conditions

- Dry ice

Product Documents

Documents for our products are available for download at takarabio.com/manuals
The following documents apply to this product:

- Yeastmaker Yeast Transformation System 2 User Manual

Quality Control Data

- Concentration measured and confirmed by OD₂₆₀.
- Size distribution measured and confirmed by agarose gel analysis.
- Functionally tested and confirmed for use in transformation with the Yeastmaker Yeast Transformation System 2 (Takara Bio, Cat. No. 630439).

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: technical_support@takarabio.com

United States/Canada
800.662.2566
(100923)

Asia Pacific
+1.650.919.7300

Europe
+33.(0)1.3904.6880

Japan
+81.(0)77.565.6999

Yeastmaker™ Carrier DNA

CATALOG NO.

630440

NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at <http://www.takarabio.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

TRADEMARKS:

©2023 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: technical_support@takarabio.com

United States/Canada

800.662.2566

Asia Pacific

+1.650.919.7300

Europe

+33.(0)1.3904.6880

Japan

+81.(0)77.565.6999

10/9/2023