pACT2 AD Vector Information

PT3022-5 Catalog 638822

GenBank Accession No.: U29899.



Restriction Map and Multiple Cloning Site (MCS) of pACT2 AD. Unique restriction sites are in bold.



United States/Canada 800.662.2566 Asia Pacific +1.650.919.7300 Europe +33.(0)1.3904.6880 Japan +81.(0)77.543.6116

Clontech Laboratories, Inc. A Takara Bio Company 1290 Terra Bella Ave. Mountain View, CA 94043 Technical Support (US) E-mail: tech@clontech.com www.clontech.com Description:

pACT2 generates a fusion of the GAL4AD (amino acids 768–881), an HA epitope tag, and a protein of interest (or protein encoded by a cDNA in a fusion library) cloned into the MCS in the correct orientation and reading frame. pACT2, which is derived from pACT (1), contains a unique *Eco*R I site in the MCS. The hybrid protein is expressed at high levels in yeast host cells from the constitutive ADH1 promoter (*P*); transcription is terminated at the ADH1 transcription termination signal (*T*). The protein is targeted to the yeast nucleus by the nuclear localization sequence from SV40 T-antigen which has been cloned into the 5' end of the GAL4 AD sequence (2). pACT2 is a shuttle vector that replicates autonomously in both *E. coli* and *S. cerevisiae* and carries the *bla* gene, which confers ampicillin resistance in *E. coli*. pACT2 also contains the *LEU2* nutritional gene that allows yeast auxotrophs to grow on limiting synthetic media. Transformants with AD/library plasmids can be selected by complementation by the *LEU2* gene by using an *E. coli* strain that carries a *leuB* mutation (e.g., HB101).

Note: The Sfi I and Sma I sites in the MCS tend to compress during sequencing.

(PR8X956)

Location of features:

- 2µ origin of replication: 1-2055
- LEU2 coding sequences Start codon (ATG): 2474–2476 Stop codon: 3566–3568
- Lox sites: Lox 1: 4268-4327; Lox 2: 4367-4412
- Transcription termination signal Fragment carrying the *S. cerevisiae* ADH1 terminator: 4415–4742 Translation stop codon: 4970–4972
- Multiple cloning site: 4927-5079
- Hemagglutinin (HA) epitope: 5042-5068
- GAL4 activation domain polypeptide Start codon (ATG): 5486–5488
 GAL4 codons 768–881: 5081–5419
 SV40 T-antigen nuclear localization signal: 5424–5478
- Promoter fragment carrying the S. cerevisiae ADH1 promoter: 5504-5901
- pBR322 plasmid replication origin: 6336-6979
- Ampicillin resistance gene

Promoter: -35 region: 8052-8057; -10 region: 8029-8034 Transcription start point: 8022 Ribosome binding site: 7995-7999 β-lactamase coding sequences: Start codon (ATG): 7985-7987 Stop codon: 7127-7129 β-lactamase signal peptide: 7919-7987 β-lactamase mature protein: 7130-7918

Primer locations:

- MATCHMAKER 5' LD-Insert Screening Amplimer: 5122–5091
- MATCHMAKER 3' LD-Insert Screening Amplimer: 4943–4974
- GAL4 AD Sequencing Primer: 5153–5137

Propagation in E. coli

- Suitable host strains: DH5 α and other general purpose strains.
- Selectable marker: plasmid confers resistance to ampicillin (50 µg/ml) to E. coli hosts.
- E. coli replication origin: pBR322
- Copy number: 15–20
- Propagation in S. cerevisiae
- Suitable host strains: Y187(α), Y190 (a), SFY526(a), CG1945(a), or HF7c(a)
- Selectable marker: LEU2
- S. cerevisiae replication origin: 2µ
- Copy number: multiple copy

Note: The attached sequence file has been compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by Clontech. This vector has not been completely sequenced.

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc. Clontech is a Takara Bio Company. ©2006