

pCherryPicker Control

Catalog No(s).

632578 & 632579 (Not sold separately). Sold as part of 632572 & 632573.

Amount

20 μ l pCherryPicker Control Vector (500 ng/ μ l)

Lot Number

Specified on product label.

pCherryPicker Control Vector Information




Figure 1. pCherryPicker Control vector map.

Location of Features

- $P_{CMV IE}$ (human cytomegalovirus immediate early promoter): 1–589
- CherryPicker (mCherry-transferrin receptor membrane-anchor domain fusion): 597–1526
- SV40 early polyA signals: 1679–1713
- f1 origin of replication: 1776–2231 (complementary)
- $P_{SV40 e}$ (SV40 early promoter): 2405–2673
- SV40 origin of replication: 2572–2710
- Kan^r/Neo^r (kanamycin/neomycin resistance gene): 2756–3550
- HSV TK polyA signal: 3786–3804
- pUC origin of replication: 4135–4778

Clontech Laboratories, Inc.

A Takara Bio Company
1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: tech@clontech.com

United States/Canada	Asia Pacific	Europe	Japan
800.662.2566	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.543.6116
Late update: July 12, 2011			

Vector Map

Cat. No(s). 632572, 632573

pCherryPicker Control

NOTICE TO PURCHASER:

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, in vitro diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with the licensing requirements listed below and described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.