



pCMV DsRed-Express2 Vector Map.

Description

pCMV DsRed-Express2 is a mammalian expression vector designed to be used for whole cell labeling or as a marker of cotransfection. In mammalian cells, the vector constitutively expresses DsRed-Express2, a variant of the *Discosoma sp.* red fluorescent protein, DsRed (1). DsRed-Express2 retains the fast maturation and high photostability characteristic of its predecessor, DsRed-Express (2), and has been engineered (through additional amino acid substitutions) for increased solubility (3). Although it most likely forms the same tetrameric structure as wild-type DsRed, DsRed-Express2 displays a greatly reduced tendency to aggregate, resulting in reduced cyto- and phototoxicity, and making DsRed-Express2 much better suited for *in vivo* applications involving sensitive cells, such as primary or stem cells. In fact, DsRed-Express2 has been shown to be the best red fluorescent protein for whole cell labeling applications (3). DsRed-Express2 also exhibits extremely low residual green fluorescence, which allows cells expressing the protein to be effectively separated from other fluorescently labeled cell populations by flow cytometry.

The DsRed-Express2 gene is positioned just downstream of the constitutively active human cytomegalovirus immediate early promoter ($P_{CMV IE}$). As a result, mammalian cells transfected with this vector will constitutively express the red fluorescent protein. A Kozak consensus sequence (4) has been placed immediately upstream of the DsRed-Express2 coding sequence to enhance translational efficiency in eukaryotic cells. SV40 polyadenylation signals downstream of the DsRed-Express2 gene direct proper processing of the 3' end of the DsRed-Express2 mRNA. The vector backbone contains an SV40 origin for replication in mammalian cells expressing the SV40 large T antigen, a pUC origin of replication for propagation in *E. coli*, and an f1 origin for single-stranded DNA production. A neomycin

(PR8Z2653; published 2 December 2008)



Clontech

United States/Canada
800.662.2566

Asia Pacific
+1.650.919.7300

Europe
+33.(0)1.3904.6880

Japan
+81.(0)77.543.6116

Clontech Laboratories, Inc.
A Takara Bio Company
1290 Terra Bella Ave.
Mountain View, CA 94043
Technical Support (US)
E-mail: tech@clontech.com
www.clontech.com

resistance cassette (Neo^r) allows stably transfected eukaryotic cells to be selected using G418 (5). This cassette consists of the SV40 early promoter, a Tn5 kanamycin/neomycin resistance gene, and herpes simplex virus thymidine kinase (HSVTK) polyadenylation signals. A bacterial promoter upstream of the cassette expresses kanamycin resistance in *E. coli*.

Use

pCMV DsRed-Express2 can be used for whole cell labeling or as a cotransfection marker. The vector can be transfected into mammalian cells using any standard transfection method. Cells expressing DsRed-Express2 (excitation and emission maxima: 554 nm and 591 nm, respectively) can be detected by fluorescence microscopy or flow cytometry 8-12 hours after transfection. If required, stable transfectants can be selected using G418. After cotransfection with pCMV DsRed-Express2 and an expression construct of interest, cells can also be sorted by flow cytometry to enrich for transfected cells.

Location of features

- P_{CMVIE} (human cytomegalovirus immediate early promoter): 1–589
- DsRed-Express2 (*Discosoma sp.* red fluorescent protein variant)
 - Kozak consensus translation initiation site: 621–631
 - Start codon (ATG): 628–630; Stop codon: 1303–1305
- Kan^r/Neo^r (kanamycin/neomycin resistance gene)
 - Neomycin phosphotransferase coding sequences: 2534–3328
- pUC origin of replication: 3913–4556

Propagation in *E. coli*

- Recommended host strain: DH5 α , HB101, and other general purpose strains. Single-stranded DNA production requires a host containing an F plasmid such as JM109 or XL1-Blue.
- Selectable marker: plasmid confers resistance to kanamycin (50 μ g/ml) in *E. coli* hosts.
- *E. coli* replication origin: pUC
- Copy number: high
- Plasmid incompatibility group: pMB1/ColE1

Excitation and emission maxima of DsRed-Express2

- Excitation maximum = 554 nm
- Emission maximum = 591 nm

References

1. Matz, M. V. *et al.* (1999) *Nat. Biotechnol.* **17**(10):969-973.
2. Bevis, B. J. & Glick, B. S. (2002) *Nat. Biotechnol.* **20**(1):83-87. Erratum in *Nat. Biotechnol.* (2002) **20**(11):1159
3. Strack, R. L. *et al.* (2008) *Nat. Methods* **5**(11):955-957.
4. Kozak, M. (1987) *Nucleic Acids Res.* **15**(20): 8125-8148
5. Gorman, C. (1985) In *DNA Cloning: A Practical Approach, Vol. II*. Ed. D. M. Glover. (IRL Press, Oxford, U.K.), pp. 143-190.

Note: The vector sequence was compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by Clontech. This vector has not been completely sequenced.

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

CMV Sequence:

DsRed-Express & DsRed-Express2:

Living Colors® Products AcGFP1, DsRed, HcRed, AsRed, AmCyan, ZsGreen, ZsYellow and their variants:

Not-For-Profit Entities: Orders may be placed in the normal manner by contacting your local representative or Clontech Customer Service at 650.919.7300. At its discretion, Clontech grants Not-For-Profit Entities a non-exclusive, personal, limited license to use this product for non-commercial life science research use only. Such license specifically excludes the right to sell or otherwise transfer this product, its components or derivatives thereof to third parties. No modifications to the protein coding sequence may be made without express written permission from Clontech. Any other use of this product requires a license from Clontech. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at licensing@clontech.com.

For-Profit Entities wishing to use this product are required to obtain a license from Clontech. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at licensing@clontech.com or [click here](#) for more information.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc., unless noted otherwise. Clontech is a Takara Bio Company. ©2008 Clontech Laboratories, Inc.