


pGAD GH Activation Domain Hybrid Vector Information

PT3104-5

GenBank Accession No.: Submission in Progress.

Cat. No. 638853



Restriction Map and Multiple Cloning Site (MCS) of pGAD GH. Unique restriction sites are in bold.

Description:

pGAD GH generates a hybrid protein that contains the sequences for the GAL4 activation domain (amino acids 768–881). pGAD GH has unique restriction sites located in the MCS region at the 3'-end of the open reading frame for the activation domain sequence. For the construction of a hybrid protein, the gene encoding the protein of interest (or a collection of cDNAs) is ligated into the MCS in the correct orientation and with the correct reading frame such that a fusion protein is generated. The fusion protein is expressed at high levels in yeast host cells from the constitutive ADH1 promoter; transcription is terminated at the ADH1 transcription termination signal. The hybrid protein is targeted to the yeast nucleus by nuclear localization sequences that have been added to the activation domain sequence from a heterologous source (Chien *et al.*, 1991). pGAD GH is a shuttle vector that replicates autonomously in both *E. coli* and *S. cerevisiae*. It carries the *bla* gene (for ampicillin resistance in *E. coli*) and the *LEU2* nutritional marker that allow yeast auxotrophs carrying pGAD GH to grow on limiting synthetic medium lacking Leu.

(PR85430)



United States/Canada
800.662.2566
Asia Pacific
+1.650.919.7300
Europe
+33.(0)1.3904.6880
Japan
+81.(0)77.543.6116

Clontech Laboratories, Inc.
A Takara Bio Company
1290 Terra Bella Ave.
Mountain View, CA 94043
Technical Support (US)
E-mail: tech@clontech.com
www.clontech.com

Location of features

- Full-length *S. cerevisiae ADH1* promoter: 6938–400
- GAL4 activation domain (AD) polypeptide with SV40 NLS: 413–985*
- Multiple cloning site: 830–903
- T7 RNA polymerase promoter/sequencing primer: 934–913
- Transcription termination signal
 - Fragment carrying the *S. cerevisiae ADH1* terminator: 1220–1547
- *LEU2* coding sequences: 2757–1666
- Col E1 plasmid replication origin: 3457–4100
- Ampicillin resistance gene: 5377–4520

* This represents the length of the polypeptide without a cloned insert. Depending on the reading frame of the cloned insert, the location of the stop codon will vary as shown in the figure on page 1.

Primer locations

- MATCHMAKER 5' AD LD-Insert Screening Amplimer (Cat No. 9103-1): 779–810
- MATCHMAKER 3' AD LD-Insert Screening Amplimer (Cat No. 9103-1): 1020–989

Propagation in *E. coli*

- Suitable host strains: DH5a, DH10 & other general purpose strains.
- Selectable marker: plasmid confers resistance to ampicillin (100 µg/ml) to *E. coli* hosts.
- *E. coli* replication origin: Col E1
- Copy number: 15–20

Propagation in *S. cerevisiae*

- Suitable host strains: Y187(α), Y190(a), SFY526(a), CG1945(a), or HF7c(a)
- Selectable marker: *LEU2*
- *S. cerevisiae* origin: 2 µ

REFERENCES

1. Bartel, P. L., Chien, C.-T., Sternglanz, R. & Fields, S. (1993) Using the two-hybrid system to detect protein-protein interactions. In *Cellular Interactions in Development: A Practical Approach*, D. A. Hartley, Ed., Oxford University Press, Oxford; pp 153–179.
2. Chien, C. T., Bartel, P. L., Sternglanz, R. & Fields, S. (1991) The two-hybrid system: A method to identify and clone genes for proteins that interact with a protein of interest. *Proc. Natl. Acad. Sci. USA* **88**:9578–9582.

Note: The attached sequence file has been compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by Clontech. This vector has not been completely sequenced.

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc.
Clontech is a Takara Bio Company. ©2006