



pHcRed1

Catalog No(s).
632410

pHcRed1 Vector Information



pHcRed1 vector map and multiple cloning site.

Clontech Laboratories, Inc.

A Takara Bio Company
1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: tech@clontech.com

United States/Canada 800.662.2566
Last update: 033012

Asia Pacific +1.650.919.7300

Europe +33.(0)1.3904.6880

Japan +81.(0)77.543.6116

Location of Features

- P_{lac} (lac promoter): 95–178
- *lacZ*-HcRed1 fusion protein expressed in *E. coli*: 217–1008
- 5' MCS (5' multiple cloning site): 234–292
- Kozak consensus sequence: 282–292
- HcRed1: 289–975
- 3' MCS (3' multiple cloning site): 992–1076
- Amp^r (ampicillin resistance gene; β -lactamase): 1522–2382
- pUC origin of replication: 2530–3173

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with the licensing requirements described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

Clontech and the Clontech logo are trademarks of Clontech Laboratories, Inc. All other marks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech is a Takara Bio Company. ©2011 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.