


pHet-Act1-2

Catalog No(s).

635082 (Not sold separately).

Sold as part of 635081.

pHet-Act1-2 Vector Information



pHet-Act1-2 vector map.

Clontech Laboratories, Inc.

A Takara Bio Company

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: tech@clontech.com

United States/Canada

800.662.2566

Last update: 071612

Asia Pacific

+1.650.919.7300

Europe

+33.(0)1.3904.6880

Japan

+81.(0)77.543.6116

Location of Features

- P_{CMV} (human cytomegalovirus promoter): 1–750
- DmrC (dimerization domain C): 1108–1392
- p65 AD (NF- κ B p65 activation domain, aa 361–551): 1399–1980
- IRES (internal ribosome entry site): 2018–2613
- NLS (nuclear localization signal): 2656–2682
- ZFHD1 (ZFHD1 DNA binding domain): 2683–3054
- DmrA (dimerization domain A): 3055–3378, 3385–3705, 3709–4032
- f1 ori (f1 origin of replication): 4410–4864
- P_{SV40e} (SV40 early promoter and enhancer): 4927–5344
- SV40 ori (SV40 origin of replication): 5243–5308
- Puro^r (puromycin resistance gene; puromycin acetyltransferase): 5388–5984
- P_{Amp}^r (ampicillin resistance promoter): 6405–6509
- Amp^r (ampicillin resistance gene; β -lactamase): 6510–7367
- pUC ori (pUC-derived origin of replication): 7541–8129

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with the licensing requirements described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

Clontech and the Clontech logo are trademarks of Clontech Laboratories, Inc. All other marks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. Clontech is a Takara Bio Company. ©2012 Clontech Laboratories, Inc.

This document has been reviewed and approved by the Clontech Quality Assurance Department.