

pMetLuc-Mem Control

Catalog No(s).
631758

Amount
20 μ l pMetLuc-Mem Control Vector (500 ng/ μ l)

Lot Number
Specified on product label.

pMetLuc-Mem Control Vector Information




Figure 1. pMetLuc-Mem Control vector map.

Location of Features

- $P_{CMV IE}$ (human cytomegalovirus immediate early promoter): 1–589
- TCA-TFR/MetLuc (chimeric non-secreted, membrane-bound *Metridia* luciferase): 597–1427
- SV40 early polyA signals: 1583–1617
- f1 origin of replication: 1680–2135 (complementary)
- SV40 origin of replication: 2476–2614
- Kan^r/Neo^r (kanamycin/neomycin resistance gene): 2660–3454
- HSV TK polyA signal : 3690–3708
- pUC origin of replication: 4039–4682

Clontech Laboratories, Inc.

A Takara Bio Company
1290 Terra Bella Avenue, Mountain View, CA 94043, USA
U.S. Technical Support: tech@clontech.com

United States/Canada 800.662.2566
Asia Pacific +1.650.919.7300
Europe +33.(0)1.3904.6880
Japan +81.(0)77.543.6116
Late update: July 12, 2011

Vector Map

Cat. No(s). 631758

pMetLuc-Mem Control

NOTICE TO PURCHASER:

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, in vitro diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Clontech Laboratories, Inc.

Your use of this product is also subject to compliance with the licensing requirements listed below and described on the product's web page at <http://www.clontech.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.