

## Mouse BCR IgG H/K/L Indexing Primer Set HT for Illumina®

Catalog Nos.	Amount	Lot Number
634428 (Not sold separately; sold as a part of 634422)	12 rxns	Specified on product label.
634429 (Not sold separately; sold as a part of 634423)	48 rxns	Specified on product label.
634430 (Not sold separately; sold as a part of 634424)	96 rxns	Specified on product label.

### Description

The Mouse BCR IgG H/K/L Indexing Primer Set HT for Illumina contains PCR primers for the amplification of indexed, paired-end Illumina-compatible sequencing libraries.

### Package Contents

<u>634428</u> (12 rxn)	<u>634429</u> (48 rxn)	<u>634430</u> (96 rxn)	
54 µl	54 µl	54 µl	BCR Primer 2V D701 (12.5 µM)
54 µl	54 µl	54 µl	BCR Primer 2V D702 (12.5 µM)
54 µl	54 µl	54 µl	BCR Primer 2V D703 (12.5 µM)
54 µl	54 µl	54 µl	BCR Primer 2V D704 (12.5 µM)
54 µl	54 µl	54 µl	BCR Primer 2V D705 (12.5 µM)
54 µl	54 µl	54 µl	BCR Primer 2V D706 (12.5 µM)
54 µl	54 µl	54 µl	BCR Primer 2V D708 (12.5 µM)
54 µl	54 µl	54 µl	BCR Primer 2V D709 (12.5 µM)
54 µl	54 µl	54 µl	BCR Primer 2V D710 (12.5 µM)
54 µl	54 µl	54 µl	BCR Primer 2V D711 (12.5 µM)
54 µl	54 µl	54 µl	BCR Primer 2V D712 (12.5 µM)
18 µl	18 µl	18 µl	mBCR Primer 2H D501 (12.5 µM)
-	18 µl	18 µl	mBCR Primer 2H D502 (12.5 µM)
-	18 µl	18 µl	mBCR Primer 2H D503 (12.5 µM)
-	18 µl	18 µl	mBCR Primer 2H D504 (12.5 µM)
-	-	18 µl	mBCR Primer 2H D505 (12.5 µM)
-	-	18 µl	mBCR Primer 2H D506 (12.5 µM)
-	-	18 µl	mBCR Primer 2H D507 (12.5 µM)
-	-	18 µl	mBCR Primer 2H D508 (12.5 µM)
18 µl	18 µl	18 µl	mBCR Primer 2K D501 (12.5 µM)
-	18 µl	18 µl	mBCR Primer 2K D502 (12.5 µM)
-	18 µl	18 µl	mBCR Primer 2K D503 (12.5 µM)
-	18 µl	18 µl	mBCR Primer 2K D504 (12.5 µM)
-	-	18 µl	mBCR Primer 2K D505 (12.5 µM)
-	-	18 µl	mBCR Primer 2K D506 (12.5 µM)

*Package Contents continued on Page 2.*

#### Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: [technical\\_support@takarabio.com](mailto:technical_support@takarabio.com)

# Certificate of Analysis

Mouse BCR IgG H/K/L Indexing Primer Set HT for Illumina  
(Not sold separately)

Cat. Nos. 634428, 634429 & 634430

Sold as a part of 634422, 634423 & 634424

---

<u>634428</u> (12 rxn)	<u>634429</u> (48 rxn)	<u>634430</u> (96 rxn)	
-	-	18 µl	mBCR Primer 2K D507 (12.5 µM)
-	-	18 µl	mBCR Primer 2K D508 (12.5 µM)
18 µl	18 µl	18 µl	mBCR Primer 2L D501 (12.5 µM)
-	18 µl	18 µl	mBCR Primer 2L D502 (12.5 µM)
-	18 µl	18 µl	mBCR Primer 2L D503 (12.5 µM)
-	18 µl	18 µl	mBCR Primer 2L D504 (12.5 µM)
-	-	18 µl	mBCR Primer 2L D505 (12.5 µM)
-	-	18 µl	mBCR Primer 2L D506 (12.5 µM)
-	-	18 µl	mBCR Primer 2L D507 (12.5 µM)
-	-	18 µl	mBCR Primer 2L D508 (12.5 µM)

## Storage Conditions

- Store at -20°C.

## Expiration Date

- Specified on product label.

## Shipping Conditions

- Dry ice

## Product Documents

Documents for our products are available for download at [takarabio.com/manuals](http://takarabio.com/manuals)

The following documents apply to this product:

- SMARTer Mouse BCR IgG H/K/L Profiling Kit User Manual

## Quality Control Data

All primers were PAGE purified and analyzed by ESI-MS.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

## Mouse BCR IgG H/K/L Indexing Primer Set HT for Illumina® - 12

### CATALOG NOS.

634428, 634429 & 634430

### NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at <http://www.takarabio.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

### TRADEMARKS:

**©2023 Takara Bio Inc. All Rights Reserved.**

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

---

#### Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: [technical\\_support@takarabio.com](mailto:technical_support@takarabio.com)

**United States/Canada**

800.662.2566

**Asia Pacific**

+1.650.919.7300

**Europe**

+33.(0)1.3904.6880

**Japan**

+81.(0)77.565.6999

3/10/2023